

Product datasheet for **RR205600L3V**

Slc25a21 (NM_133614) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc25a21 (NM_133614) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc25a21
Synonyms:	Odc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133614
ORF Size:	894 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205600).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133614.2 , NP_598298.1
RefSeq Size:	2470 bp
RefSeq ORF:	897 bp
Locus ID:	171151
UniProt ID:	Q99JD3
Cytogenetics:	6q23
Gene Summary:	mitochondrial oxodicarboxylate carrier; involved in the transport C5-C7 oxodicarboxylates across the inner membranes of mitochondria [RGD, Feb 2006]



[View online »](#)