

Protein Sequence: >RR205607 representing NM_001109139
 Red=Cloning site Green=Tags(s)


MASATAAAAPGEAEETTRLRKPRFSFEENQILIREVRAHYPQLYGAQSRVSVAEERRRVWDSIATKINGI
 TSWKRTGQEVQKRWNDFKRRTKEKLARVPHSTQGAGPAEDAFSAEETIFAILGPGVAGPGAGAGAEES
 PAAASSQPQASTASTQRYVLSDDRQDRRADTPAQSKGSSSPESWARPSCNPQEAKERESPSPAA MQPV
 QLPRLALSPPLPAPPPPTALAQVAPSSPTPPRPTSAPEQSLDFLRAQQETANAIRELAGTLRQGLAK
 LSEALSALLPLPGTPADPQPPPPPPPPPPPKPVLPPPAPKVELAPEVSVVAAVVDGAVVAARGVIIS
 PRSEEGVPKPPPVAPLPLHDSPPHKRRKGFPTRKRRGRWKSP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001109139

ORF Size: 1176 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001109139.1, NP_001102609.1

RefSeq Size: 1866 bp


RefSeq ORF: 1179 bp

Locus ID: 499090

Cytogenetics: 1q21

MW: 41.7 kDa

Product images:



Circular map for RR205607