

Product datasheet for **RR205612**

Tmeff2 (NM_001108795) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmeff2 (NM_001108795) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tmeff2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR205612 representing NM_001108795
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGCTGTGGGAGTCCCCGCGGCAGTGCAGCAGCTGGACACTTTGCGAGGGCTTCTGCTGGCTGCTGC
TGCTGCCGGTGACGTGCTTATCATAGCCCGCCGGTGAAACTCGCTGCTTCCCTACCTCCTTAAGTGA
CTGCCAAACGCCACCGGCTGGAAGTCTCCGGTTACGATGACAGAGAAAATGATCTCTTCTGTGTGAT
ACCAACACCTGTAATTTGACGGGAATGTCTAAGGATTGGAGACTGTGACTTGCCTGTGTCAGTTCA
AGTGTAACAGTGACTACGTGCCTGTGTGGCTCCAATGGAGAGAGCTACCAGAATGAATGCTACCTGAG
GCAAGTGCCTGCAAGCAACAGAGTGAGATACTTGTGGTCTCAGAAGGATCCTGTGCTACAGATGCAGGA
TCGGGATCTGGAGACGGAGTCCATGAAGGCTCTGGAGAGACCAGTCAGAAGGAGACATCTACCTGTGACA
TTTGCCAGTTTGGTGCGGAATGCGATGAAGATGCTGAGGATGTCTGGTGTGTGCAACATTGACTGTTT
TCAAACCAACTTAAATCCCCTCTGTGCTTCTGATGGGAAATCTTATGATAATGCATGTCAAATCAAAGAA
GCGTCATGTGAGAAACAGGAGAAAATTGAAGTCATGTCTCTGGGTCGATGTCAAGATAACACAACACTCAA
CCACTAAATCTGAAGATGGGCATTATGCAAGAACAGATTATGCAGATGTGATGCTGGTTACTGGTCAA
CACTGTGAAAAGAAGGACTACAGTGTTCTCTACGTTGTTCTGGTCTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR205612 representing NM_001108795
Red=Cloning site Green=Tags(s)

MVLWESPRQCSSWTLCEGFCWLLLLPVTLLIARPVKLAAPFPTSLSDCQPTPTGWNCSGYDDRENDLFLCD
 TNTCKFDGECLRIGDVTVCVQFKCNSDYVPVCGSNGESYQNECYLRQAACKQQSEILVSEGCATDAG
 SSGSDGVHEGSGETSQKETSTCDICQFGAECDEDAEDVWCVCNIDCSQTNFNPLCASDGTKSYDNACQIKE
 ASCQKQEKIEVMSLGRCDQNTTTTTKSEDGHYARTDYADVMLVILVNTVKRRRTTVFSTLFLVL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108795

ORF Size: 819 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001108795.1 , NP_001102265.1
RefSeq Size:	3252 bp
RefSeq ORF:	822 bp
Locus ID:	363228
Cytogenetics:	9q22
MW:	30 kDa