

Product datasheet for **RR205613L3V**

Pmfbp1 (NM_134393) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pmfbp1 (NM_134393) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pmfbp1
Synonyms:	ODF3; Stap
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_134393
ORF Size:	2913 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205613).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134393.2 , NP_599220.1
RefSeq Size:	3307 bp
RefSeq ORF:	2916 bp
Locus ID:	171414
UniProt ID:	Q9Z221
Cytogenetics:	19q12
Gene Summary:	coiled-coil protein that is a component of the outer dense fibers of the sperm tail; may play a role in cytoskeleton organization [RGD, Feb 2006]


[View online »](#)