

Product datasheet for **RR205632**

Lrrc10b (NM_001107577) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc10b (NM_001107577) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lrrc10b
Synonyms:	RGD1564983
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR205632 representing NM_001107577 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCATCGCAGAGTCCACGCCGACGAGCTGCCGTCGGACGCCGAGGAGCAGCTACGCAGCGCGAAC
AGCAGCTGGAGCTGAGCGGGCGGCGTTGCGCGGCTGCCAGCGCGTGTGCGCACTGAGCCGCTTGCA
GAAGCTGTATGTAGCGGCACCGGGCTGCGCGAGCTTCCGGAGGAGATAGAGGAGCTGCGCGAGCTGCGC
ATCCTGGCGCTCGACTTCAACAAGCTGGAGCGCCTGCCGACGGCCTGTGTCGCCTGCCGCGCTCACGC
GCCTCTACCTGGGCGGAACCGTCTGCTGGCGCTGCCATCCGACTTTGCGCAGCTGCAGAGCCTGCGCTG
CCTCTGGATCGAAGGCAACTTCTGCGGCGCTTTCCGCGGCTCTGCTGCGCCTGGTGGCGCTGCAATCG
CTGCAGATGGGCGACAACCGGCTGCGCGCGCTGCCGCGGAGCTGCCGCGCATGACGGGCTTGCCTGGCC
TCTGGCTCTATGGCAATCGCTTCAAGAGTTCCCGCCCGCGCTGTTGCGCATGGGCCGGCTGCACATCCT
CGACCTCGACCGCAACCGCTGGGCGGCTTCCCGACCTGCACCCTCTGCGTGCGCTCCGAGTTTCTCC
TACGACCACAACCCAGTCACTGGACCCCGCGGTGGCAGACACCGTTTTCTAGTGGGTGAGGGCGCCG
TCGAGCGCATGGCTGAGCGTGACGAACCATTCGCGGCCGCCAGCCCGCGTCCAGCCCGGCCCTTTGA
TGATGAGGAGGAAGAAGACTTGCTCATAGGAAGTGCAGGACCCAGAGCCCTGGGGCCGCCAGGGATAAC
CTCCGAGCCTTGAAGCCCTCCAGGATTGGGTACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »


Protein Sequence: >RR205632 representing NM_001107577
 Red=Cloning site Green=Tags(s)

MGIAESTPDELPSDAEEQLRSGEQQLLESGRRLRLPSAVCALSRQLKYVSGTGLRELPEEIEELREL
 ILALDFNKLRLPDGLCRLPRLTRLYLGNRLLALPSDFAQLQSLRCLWIEGNFLRRFPRPLRLVALQS
 LQMGNRLRALPAELPRMTGLRGLWLYGNRFEEFPPALLRMGRLHILDLDNRNLGGFPLHPLRALRVFS
 YDHNPTVGPPRVADTVFLVGEGAVERMAERDEIPRPPARRPARAFDDEEEEDLLIGSAGPRALGPARDN
 LRALEAPPGLGT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001107577

ORF Size: 876 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001107577.1, NP_001101047.1

RefSeq Size: 1229 bp


RefSeq ORF: 879 bp

Locus ID: 309208

Cytogenetics: 1q43

MW: 32.9 kDa

Product images:



Circular map for RR205632