

Product datasheet for **RR205633L3V**

Got1 (NM_012571) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Got1 (NM_012571) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Got1
Synonyms:	ASAT; Aspat; cAspAT; cCAT; Gaspat
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012571
ORF Size:	1239 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205633).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012571.2 , NP_036703.2
RefSeq Size:	1998 bp
RefSeq ORF:	1242 bp
Locus ID:	24401
UniProt ID:	P13221
Cytogenetics:	1q54
Gene Summary:	cytosolic aspartate aminotransferase isozyme; enzyme activity in liver and kidney increase in response to glucocorticoids [RGD, Feb 2006]



[View online »](#)