

Product datasheet for **RR205634**

LOC498971 (NM_001047936) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC498971 (NM_001047936) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LOC498971
Synonyms:	Rtap2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR205634 representing NM_001047936 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAATCAGGCAATGCATCCATCCACCATTAGTAAGAAATTGACCAGTGTTCCAGAATTACCTTTTAAGA
AAGGACTACTGAATTTATCACCTGCACTGAAAGAAAAACATAGTCTCAAGTCTACAGATGTTAAATTTGA
ACCAATGGTCCTCAGGTCTCCCCACAGGAGAGTCGATTGTACGGTTTGCCTTGCCAATACCACTGAGA
AAGACGAAGGAGTTAATATCAAGTCATGAAATGGTCAGGAAGATTACCCAAAACCTGAAGATGGTTGTTT
CTACCTTAGAGAAGACCTTTGGAAATATCAATGAAGATAAAGAAAAAGAGCAGCAAAACCTGAAATAGA
AGGCTTCAGACAAGAAGGCTTGTCTGTTGGTGATGATGTGAAGTCATTCTTGTGTGCTGTTTCAGAAATTT
ACAGCTCAGCTGGAGGAAGCAGTTAAGGAAGAAGCTGGTGCCTTGAATCTCTTTACAAGTGGTTTCAGC
AACAAAGTCAACCAACTGGAGGAGATCGGCAAGATCAAAGTAATTTAGAAGCAGAAGTTCATCAGATAG
AAAATCAATTTTGAACATTGTACAAATTGCAAACTGGCACGTAATTTGAAGATATTGAAGGCCGCCTA
AGAGAGCGCAAGGCGCCATGCTGGCTAAAGAAGAGGATAAAGCAACTGTATCAGAATCATTGAAAAATT
ATCAACTTATAGAAAAGCAAATTGAAGAGTTCATAAAATCCCATTGATATAGAATTTCAAATTTGTTTC
TGAACTGAGTCAGGTTTGGAAACAGCCAAGTGGCTCAACAGGAACTCCTTCAGTGACTAATCGGATGAAT
ACAATGATGAAAAATTTGAAAAACAGCAACTATGTTAGAAAAAGCTTTAAGTGATCAGACTGAAGCTG
AGGCTAAATATAAAAAAATGGAACATAATGTTCAACTGTTAGTAATGGAGAAAGCCCTCCTGAAAGGTGA
AATACGAAGGCTGAGAGAAATACAGAAAGCAAAGGCTACAAGCAAAGAAGTGCGACCCAAGAAATCTGGA
AAATCAGAAAAGAAAAAGATAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RR205634 representing NM_001047936
 Red=Cloning site Green=Tags(s)


MNQAMHPSTISKKLTSPPELPFKKGLLNLSPALKEKHSKSTDVKIEPMVLRSPPTGESIVRFALPIPLR
 KTKELISSHEMVRKITQNLKMOVSTLEKTFGNINEDKEKRAAKPEIEGFRQEGLSVGDDVKSFLCCSEF
 TAQLEEAVKEERGALESYKWFQQQVNLLEEIGKQSNLEAEVQSDRKSILNIVQIAKLARKFEDIEGRL
 RERKAAMLAKEDKATVSESLKNYQLIEKQIEEFIKSHSDIEFQIVSETESGLEQPSGSTGTPSVNRMN
 TMMKIFENQTTMLEKALSDQTEAEAKYKKMETNVQLLVMEKALLKGEIRRLREIQKAKATSKEVRPKKSG
 KSEKKKDR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001047936

ORF Size: 1074 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001047936.1, NP_001041401.1


RefSeq Size: 1649 bp

RefSeq ORF: 1077 bp

Locus ID: 498971

MW: 40.7 kDa

Product images:



Circular map for RR205634