

Product datasheet for **RR205665L3V**

Jdp2 (NM_053894) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Jdp2 (NM_053894) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Jdp2
Synonyms:	Jundm2; Jundp2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053894
ORF Size:	489 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205665).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053894.1 , NP_446346.1
RefSeq Size:	878 bp
RefSeq ORF:	492 bp
Locus ID:	116674
UniProt ID:	Q78E65
Cytogenetics:	6q31
Gene Summary:	transcription factor that represses AP-1-mediated activation [RGD, Feb 2006]



[View online »](#)