

## Product datasheet for **RR205694**

### **Tnip1 (NM\_001108826) Rat Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Tnip1 (NM_001108826) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tnip1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>RR205694 representing NM\_001108826  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGAAGGGAGAGACCCTACCGGATCTACGACCCAGGGGCGAGCAGCCTCTGGGAGAGGAGGTGTCCG  
CAGCTTTTGACCGTCTAGTGGAGGAGAATACTCGGCTGAAGGGAAAAATACAAGGGATAAAGATGTTAGG  
AGAGCTTCTGGAGGAATCTCAAATGGAAGCGTCCAGACTCCGGCAGAAGGCGGAGGCGCTGGTCAGGGAC  
AGCGAGCTGCCTCCACCACCACCTGCTCCCTCCTTGGTCTCCTTTGATCCCCTGGCTGAGCTCACAGGAC  
AGGATACAAAGGTTCAAGGCACATCCTACTACTACCACCACCACCACCACCACCACCACCACCACCACCAT  
TGGAACTCCATGGAGAAGTCAAGTGTGAGCTCCAGCTCCATCCAATGGCGCTCCTCTGACTTTGAA  
GTGGTCCCTCCTGAGGAGCAGAGTTCACCTGAACTGGCCACCCTCTAGGAACGTGATGGAGCTGGGGT  
CCCCACCCAGAGGACGCAATCTGATGCTCCACCTGCAGCGCTGGAGACCACCCTGAGTGTGTGTGC  
AGAGGAACCAGATCACAGCCAGCTTTCAACCCTGGGCGCATGGCCCTCGAGTTCAACAGGTTGGCC  
TCCAAAGTGCATAAAAAATGAGCAACGTACTCCATCCTGCAGACCTTATGTGAGCAGCTGCGCCAGGAGA  
ATGAAGCCCTGAAGGCCAAGCTGGACAAGGCCTGGAACAGCGGGATCAGGCTGTGAGAGGCTGCGGGGA  
GGAAAAACAGGAGCTCAAGAGACTGTTGATGAACAGCAACTGCAAGAGGGACTCTGTGGGAGCCAGC  
TCCCCAAAGACAGAGGGTGTGGCAAGAAGGCGTGGCTGGACAGCAGCAGGCCAGTATGACGGCAAGTA  
AGGTCCCTGAAGCAGGGGCTTCGGTACAGCTGAGAAGAAGGTGAAGTTGCTGGAACAGCAGCGCATGGA  
GCTGCTGGAAGTGAACAAGCAGTGGGACCAGCATTTCCGGACATGAAGCAACAGTATGAGCAGAAGATC  
ACAGAGCTTCGCCAGAAGCTGGTGGACCTGCAGAAACAGGTGACTGAGCTGGAGGCTGAGCGGAGCAGA  
AGCAGCGTGACTTTGACCGGAACTCCTCCTGGCTAAATCTAAGATTGAGATGGAAGAGACTGACAAGGA  
ACAGCTGACAACAGAGGCCAAGGAACTTCGCCAGAAGTCAAGTACCTGCAGGATCAGCTGAGCCACTT  
ATGAGGCAACGAGAATACCAGGAAAAGGAGATCCAGCGGCTCAACAAGGCCCTGGAGGAGCCCTCAGCA  
TCCAGGCTCTCCCTCATCCCCACTGCAGCTTTTGGGAGCCAGAAGGGTTCGGGGCCATCTGAGGAA  
GCAGGAACTGGTACGCAGAATGAGCTATTGAAACAGCAGGTAAGATCTTTGAAGAGGACTTCAGAGG  
GAACGGAGTGACCGTGAAGCATGAACGAAGAGAAGGAGGAGCTGAAGAAGCAGGTGGAGAAGCTGCAGG  
CCCAGGTACCCCTGACCACTGCCAGCTCAGAACACTCAAAGAGGAGGAGAAGGCCAAGGAAGCCCTCAA  
ACAACAGAAGAGGAAGGCAAAGGCTCAGGAGAGCGCTACCATGTGGAACCCACCCTGAGCACCTTTGT  
GGCGCTATCCCTATGCCTACCCACCCATGCCAGCCATGGTACCTCACCATGCATACAAGGACTGGTCCC  
AGATCCGATACCCCTCACCCCTGTGCCATGGAGCACCCGCCCCACACCCCAACTCTCGTCTCTTCCA  
TCTGCCAGAGTACACCTGGCGTCCACCCTGTGCAGGAATTCGGAATCAGAACTCTCAAGTATGACCCG  
CCCCAGACAGGCTGCAGAGCCAGAGTCTGCAGACAATGACTGTGACGGGCCCCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTAA

**Protein Sequence:** >RR205694 representing NM\_001108826  
 Red=Cloning site Green=Tags(s)

```
MEGRGPYRIYDPGGSTPLGEEVSAADFRLVEENTRLKGIQGIKMLGELLEESQMEASRLRQKAEALVRD
SELP PPPPPAPSLVSFDPLAELTGQDTKVQAHPTTTTTTTTTTTTTTTIGNSMEKSVSASKPPSNGASSDFE
VVPPEEQSSPETGHHSRNVMELGSPPPEDGNLMLHLQRLETTLSVCAEEDHSQFLTHLGRMALEFNRLA
SKVHKNEQRTSILQTLCEQLRQENEALKAKLDKSLEQRDQAAERLREENTELKRLLMNSNCKEGLCGQPS
SPKTEGAGKKGAVAGQQQASMTASKVPEAGAFGTAEKVKLLEQQRMELLEVNKQWDQHFRTMKQQYEQKI
TEL RQKLVDLQKQVTELEAEREQKQRFDRKLLLAKSKIEMEETDKEQLTTEAKELRQKVLYLQDQLSPL
MRQREYQEKEIQRLNKALEEALS IQASPSPAAFGSPEGVGGHLRKQELVTQNELLKQVYKIFEEEDFOR
ERSDRERMNEEKEELKKQVEKLQAQVTL TTAQLRTLKEEEKAKEALKQQRKAKASGERYHVEPHPEHLC
GAYPYAYPPMPAMVPHHAYKDWISQIRYP PPPVPMEHPPHPNSRLFHLPEYWRPPCAGIRNQNSQVMDP
PPDRPAEPESADNDCDGPQ
```


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_001108826

**ORF Size:** 1947 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM\_001108826.1, NP\_001102296.1  
RefSeq Size: 2867 bp  
RefSeq ORF: 1950 bp  
Locus ID: 363599  
Cytogenetics: 10q22  
MW: 73.5 kDa