

Product datasheet for **RR205746L3V**

Azi2 (NM_001025705) Rat Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Azi2 (NM_001025705) Rat Tagged ORF Clone Lentiviral Particle |
| Symbol: | Azi2 |
| Synonyms: | MGC112644 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001025705 |
| ORF Size: | 1212 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RR205746). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001025705.1 , NP_001020876.1 |
| RefSeq Size: | 2017 bp |
| RefSeq ORF: | 1215 bp |
| Locus ID: | 316051 |
| UniProt ID: | Q4KMA0 |
| Cytogenetics: | 8q32 |



[View online »](#)

Gene Summary:

Adapter protein which binds TBK1 and IKBKE playing a role in antiviral innate immunity. Activates serine/threonine-protein kinase TBK1 and facilitates its oligomerization. Enhances the phosphorylation of NF-kappa-B p65 subunit RELA by TBK1. Promotes TBK1-induced as well as TNF-alpha or PMA-induced activation of NF-kappa-B. Participates in IFNB promoter activation via TICAM1 (By similarity).[UniProtKB/Swiss-Prot Function]