

Product datasheet for **RR205768L4V**

Tmed2 (NM_031722) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tmed2 (NM_031722) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tmed2
Synonyms:	Rnp24
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031722
ORF Size:	603 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205768).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031722.2 , NP_113910.1
RefSeq Size:	1988 bp
RefSeq ORF:	606 bp
Locus ID:	65165
UniProt ID:	Q63524
Cytogenetics:	12q15
Gene Summary:	may play a role in intracellular vesicle-mediated protein transport from the endoplasmic reticulum to the Golgi apparatus [RGD, Feb 2006]



[View online »](#)