

Product datasheet for **RR205779L3V**

Btg2 (NM_017259) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Btg2 (NM_017259) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Btg2
Synonyms:	Agl; An; an-1; Pc3; Tis21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017259
ORF Size:	474 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205779).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017259.1
RefSeq Size:	2519 bp
RefSeq ORF:	477 bp
Locus ID:	29619
UniProt ID:	P27049
Cytogenetics:	13q13
Gene Summary:	an anti-proliferative protein; interacts with Pick1 and may have a role in PKC-mediated extracellular signal transduction and cellular differentiation [RGD, Feb 2006]



[View online »](#)