

Product datasheet for **RR205809L4V**

Ube2n (NM_053928) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2n (NM_053928) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2n
Synonyms:	MGC93937
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053928
ORF Size:	456 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205809).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053928.2 , NP_446380.1
RefSeq Size:	1575 bp
RefSeq ORF:	459 bp
Locus ID:	116725
UniProt ID:	Q9EQX9
Cytogenetics:	7q13
Gene Summary:	mediates ubiquitin-dependent proteolytic pathway; forms a thiolester linkage between the active site cysteine and the ubiquitin, supports conjunction with a third protein, ubiquitin-protein ligase (E3) [RGD, Feb 2006]



[View online »](#)