

Product datasheet for **RR205813L4V**

Ces1d (NM_133295) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ces1d (NM_133295) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ces1d
Synonyms:	Ces3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133295
ORF Size:	1695 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205813).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133295.3 , NP_579829.3
RefSeq Size:	1938 bp
RefSeq ORF:	1698 bp
Locus ID:	113902
UniProt ID:	P16303
Cytogenetics:	19p11
Gene Summary:	enzyme that hydrolyzes long-chain fatty acids and thioesters; may play a role in lipid metabolism and/or the elimination of toxic substances [RGD, Feb 2006]



[View online »](#)