

Product datasheet for **RR205822L4V**

Adipoq (NM_144744) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adipoq (NM_144744) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adipoq
Synonyms:	Acdc; Acrp30
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_144744
ORF Size:	732 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205822).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144744.3 , NP_653345.1
RefSeq Size:	1487 bp
RefSeq ORF:	735 bp
Locus ID:	246253
Cytogenetics:	11q23
Gene Summary:	adipocyte complement-related factor that plays a role in obesity-related insulin resistance [RGD, Feb 2006]


[View online »](#)