

## Product datasheet for **RR205825**

### Calhm1 (NM\_001109168) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Calhm1 (NM\_001109168) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Calhm1  
**Synonyms:** Fam26c; RGD1564891  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR205825 representing NM\_001109168  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGACAAGTTTCGGATGATCTTCCAGTCTTGCAGTCCAACCAAGAGTCCTTCATGAACGGCATCTGCG  
GCATCATGGCGCTGGCCAGTGCCAGATGACTCTGCCTTGACTTCAACTGCCCTGCTTGCCCGGCTA  
CAATGTGGTCTACAGCATGGGCATACTGCTGACGCCTCCCTGGTCTTCTCTGCTTGGCTGGTCATG  
AACAAACAATATCCATGCTAGCTGAAGAGTGAAACGCCCCCTGGGGCGCCGGGCAAGGACCCCTGCCG  
TGCTACGCTACATGTTCTGTTCCATGGCCAGCGAGCGCTCATCGCCCTGTGGTCTGGTGGCAGTCAC  
ACTGCTGGACGGCAAGTGTTCCTCTGTGCTTTCTGCACAGCCGTGCCGGTGGCCACGCTAGGCAATGGG  
AGCCTGGTGGCCGGCCTGCCTGCTCCAGAACTTGTCTGCCTGCTGGCTCGTGACCTTGCCTGAAATCT  
ACGATGGAACTGGCTGCTGGCCGAGAGGTGGCCGTGCGCTATTTGCGCTGCATCTCAGGCGCTGGG  
CTGGTCTTCTCAAGAGCAAGTACTGGTCCCACTACATTGACATTGAGCGTAAGCTTTCGATGAGACT  
GCACGGAGCATGCCAAAGCCTTTGCTAAGGTATGTATCCAGCAGTTCTTCGAAGCCATGAACCATGACCT  
GGAGCTGGGTATACCCACGGAGTACTGGCCACAGCTTCGGCCAGAACTACAGCCACAGAGGCCATCCAA  
AGTCCCTCCGACAGGACCGAAGAAGAGAGGGAGAAGCTGCGAGGCATCACTGACCAAGGCAGCATGAACA  
GGCTACTCACAAGCTGGCACAAATGCAAACCACCTTGAGGCTGGCCAGGAGGCACCACTGATGAGCAA  
TGGCTGGGCTGGGAGCGAGCCCCGGCCTCCACGCAAGGAAGTGGCCACCTACTTCAGCAAAGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR205825 representing NM\_001109168  
Red=Cloning site Green=Tags(s)

MDKFRMIFQFLQSNQESFMNGICGIMALASAQMYSAFDNCPCLPGYNVVYSMGILLTPPLVLFLLGLVM  
 NNNISMLAEWKRPLGRRAKDPAVLRMFCSMAQRALAPVVVAVTLLDGKCFLLCAFCTAVPVATLGNQ  
 SLVPGLPAPELARLLARVPCPEIYDGNWLLAREVAVRYLRCSQALGWSFVLLTLLAFVVRVSRPCFTQ  
 VAFLLKSKYWYSHYIDIERKLFDETCTEHAKAFKVCIQQFFEAMNHDLELGHHTGVLATASARTTATEAIQ  
 SPSDRTEEREKLRGITDQGSMNRLTSHWKCKPPLRLGQEAPLMSNGWAGSEPRPRPRKEVATYFSKV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001109168

**ORF Size:** 1044 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001109168.1</a> , <a href="#">NP_001102638.1</a>
<b>RefSeq Size:</b>	1047 bp
<b>RefSeq ORF:</b>	1047 bp
<b>Locus ID:</b>	499367
<b>Cytogenetics:</b>	1q54
<b>MW:</b>	39 kDa