

Product datasheet for **RR205828**

Ncbp1 (NM_001014785) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ncbp1 (NM_001014785) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ncbp1
Synonyms:	MGC105807
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RR205828 representing NM_001014785
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCGGGAGGCGGCACAGCTACGAGAACGATGGTGGACAACCTCACAAAAGGAGGAAGACATCAGATG
 CAAATGAAACGGAAGATCATTGGAACTTTAATATGCAAAGTAGGGGAAAAGAGCGCTTGCTCTTTAGA
 GAGCAATCTCGAAGGCCTGGCTGGTGTTTAGAAAGCTGATCTTCTAACTACAAGAGCAAGATCTTAAGG
 CTCCTTTGTACAGTTGCACGTCTGTACCTGAGAAGCTAAACAATTTATACAACATTAGTTGGACTGCTGA
 ATGCCAGAAATTATAACTTTGGAGGAGAATTTGTAGAAGCTATGATTCGCAACTAAAGAATCGCTAAA
 AGCAAACAATAATGAAGCTGTATTTGGTCCGTTTTCTGTCTGATCTTGTGAAGTGCATGTGATT
 GCTGCCCATCTATGGTTGCTATGTTTGAGAATTTGTAAAGTGAAGTCAAGGAAGAAGATGTGCCTCAAG
 TGGCAGCGGATGGTATGTGTACGCTTTCTGTATCCTTCCCTGGGTTGGGAAGGAGTTGATGAAAA
 GAAAGATGCAGAAATGGACCGCATCTTCTCCAGCACCGAAAGCTATCTCAAAGACGCCAGAAGACTCAT
 GTGCCATGTTACAAGTATGGACTGCCGATAAGCCACACCCCCAAGAAGATTTAGATTGCCTGTGGG
 CCCAGATTGAGAAATGAAAAAGGACCGGTGGCAGGAAAGGCACATCCTAAGACCTACCTTGCCTTTGA
 CAGCATCCTCTGTGAAGCACTACAGCACAACCTGCCTCCTTTCACACCACCGCCTCACACCGAAGACTCC
 GTGTACCAATGCCAGGGTCATCTTCAGGATGTTTACTACAGACAGCCAGAGGGTCTGTATGC
 CAGGGAGCCATTAGTAGAGAGGTTTGTATAGAGGAGAATCTTCATTGCATCATAAAATCCTACTGGAA
 GGAAAGGAAAACCTTGCTGCCAGTTGGTGTATCCAGGCAAGAACAAGATCCATTGAACTACCAC
 ATAGTGGAGGTGATCTTGCAGAGCTGTTTCCAGCTTCCAGCGCCTCCGCACATTGATGTGATACACGA
 CACTGCTGATCGAAGTGTGCAAACTTCAGCCGGCTCTCTGCCCAAGTGTGTCACAGGCCACTGAGAT
 GCTGTACATGCGTCTGGACACCATGAGTACCACGTGTGTTGACAGGTTTATTAATTGGTTTTCTCATCAT
 CTAAGTAAATTTCAATCCGCTGGAGCTGGGAAGATTGGTGTGATTGCTTACTCAGGATCTTGAAGTGC
 CCAAACCAAAGTTGTAAGAGAAGTTCTAGAGAAATGTATGAGGTTGTCTTACCATCAGCATATATTAGA
 TATCGTTCCTCCACCTTCTCAGCTCTGTGCTGCAAACCAACCTGCATTTACAAGTATGGAGACGAA
 AGTAGCAATCTTCTGGACATTCGGTGGCACTCTGTTTATCTGTTGCTTTTAAAAGTAAAGGCAACCA
 ACGATGAAATCTCAGCATTCTGAAGGACGTACCAAATCCTAACCAGGTTGATGATGACGATGAAGGGTT
 CAGATTTAATCCCTTGAAGTCTTGTCCAGACTCTGCTGCACTTAGCTGCAAGTCTTTCAGT
 CACTCCTTCAGTCTTGTCAAAGTTTTCATGAAGTCTTCAAACCTCTAGCAGAAAGTATGAAAGGGAAGC
 TGATGTGCTGAGAGTCATGTTTGAAGTCTGGAGGAATCACCCGAGATGATTGCTGTACTGGTTGATAA
 GATGATTTCGTACACAGATTGTTGACTGTGCTGCAGTAGCAAACCTGGATCTTCTCATCAGAGCTGTCTCGT
 GACTTTACGAGGTTGTTTGTGTGGGAAATTTGCACTCTACAATTCGTAAGATGAATAAGCATGTTCTGA
 AGATCCAGAAAGAGCTGGAAGAAGCTAAAGAGAAACTAGCAAGGCAGCACAAACGACGCAGTGTATGA
 TGACCGGGGCAGCGACAGGAAGGACGGTGCCTGGAGGAGCAAATCGAAAGGCTGCAGGAGAAAGTGGAG
 TCTGCTCAGAGCGAGCAAAGAACCTTTTCTCGTCATCTTCCAGCGTTTCATCATGATCTTGACTGAAC
 ACTTGGTGCCTGTGAAACCGATGGGACCAGTATATTGACCCCGTGGTATAAGAACTGCATCGAGAGGCT
 GCAGCAGATCTTCTACAGCATCACCAGATAATCCAGCAGTACATGGTGACCCTGGAGAAACCTGCTCTTC
 ACGGCAGAGTTAGATCCTCACATCCTGGCTGTGTTCCAGCAGTCTGTGCTCTGCAGGCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RR205828 representing NM_001014785
Red=Cloning site Green=Tags(s)

MSRRRHSYENDGGQPHKRRKTSANETEDHLES LICKVGEKSACSLASNLEGLAGVLEADLPNYKSKILR
LLCTVARLLPEKLTIIYTTLVGLLNARNYNFGGEFVEAMIRQLKESLKANNYNEAVYLVRFLSDLVNCHVI
AAPSMVAMFENFVSVTQEEDVPQVRRDYYVYAF LSSLPWVGKELYEKKAEMDRIFSSTESYLKRRQKTH
VPMLQVWTADKPHPQEEYLDCLWAQIQK LKDRWQERHILRPYLAFDSILCEALQHNLPPFTPPPHTEDS
VYPMPRVIFRMFDYTDDEPGVPMPGSHSVERFVIEENLHCIIKSYWKERKTCAAQLVSYPGKNKIPLNYH
IVEVIFAELFQLPAPPHIDVMTLLIELCKLQPGSLPQVLAQAATEMLYMRLDTMSTTCVDRFINWFSHH
LSNFQFRWSWEDWSDCLTQDLESPKPKFVREVLEKCMRLSYHQHILDIVPPTFSALCPANPTCIYKYGDE
SSNSLPGHSVALCLSVAFKSKATNDEIF SILKDVNPNQVDDDEGFRFNPLKIEVFVQTLHLAAKSFS
HSFSALAKFHEVFKTLAESDKGKLHVL RVMFEVWRNHPQMI AVLVDKMIRTQIVDCAAVANWIFSSLSR
DFTRLFVWEILHSTIRKMKNHVLKIQKELEEAKELARQHRRSDDDRGSDRKGAL EEQIERLQEKVE
SAQSEQKNLFLVIFQRFIMILTEHLVRCETDGT SILTPWYKNCIERLQQIFLQHHQIIQQYMTLENLLF
TAELDPHILAVFQQFCALQA

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001014785

ORF Size: 2370 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001014785.1](#), [NP_001014785.1](#)

RefSeq Size: 2850 bp

RefSeq ORF: 2373 bp

Locus ID: 298075

UniProt ID: [Q56A27](#)


Cytogenetics: 5q22

MW: 91.9 kDa

Gene Summary:

Component of the cap-binding complex (CBC), which binds cotranscriptionally to the 5'-cap of pre-mRNAs and is involved in various processes such as pre-mRNA splicing, translation regulation, nonsense-mediated mRNA decay, RNA-mediated gene silencing (RNAi) by microRNAs (miRNAs) and mRNA export. The CBC complex is involved in mRNA export from the nucleus via its interaction with ALYREF/THOC4/ALY, leading to the recruitment of the mRNA export machinery to the 5'-end of mRNA and to mRNA export in a 5' to 3' direction through the nuclear pore. The CBC complex is also involved in mediating U snRNA and intronless mRNAs export from the nucleus. The CBC complex is essential for a pioneer round of mRNA translation, before steady state translation when the CBC complex is replaced by cytoplasmic cap-binding protein eIF4E. The pioneer round of mRNA translation mediated by the CBC complex plays a central role in nonsense-mediated mRNA decay (NMD), NMD only taking place in mRNAs bound to the CBC complex, but not on eIF4E-bound mRNAs. The CBC complex enhances NMD in mRNAs containing at least one exon-junction complex (EJC) via its interaction with UPF1, promoting the interaction between UPF1 and UPF2. The CBC complex is also involved in 'failsafe' NMD, which is independent of the EJC complex, while it does not participate in Staufen-mediated mRNA decay (SMD). During cell proliferation, the CBC complex is also involved in microRNAs (miRNAs) biogenesis via its interaction with SRRT/ARS2 and is required for miRNA-mediated RNA interference. The CBC complex also acts as a negative regulator of PARN, thereby acting as an inhibitor of mRNA deadenylation. In the CBC complex, NCBP1/CBP80 does not bind directly capped RNAs (m7GpppG-capped RNA) but is required to stabilize the movement of the N-terminal loop of NCBP2/CBP20 and lock the CBC into a high affinity cap-binding state with the cap structure. Associates with NCBP3 to form an alternative cap-binding complex (CBC) which plays a key role in mRNA export and is particularly important in cellular stress situations such as virus infections. The conventional CBC with NCBP2 binds both small nuclear RNA (snRNA) and messenger (mRNA) and is involved in their export from the nucleus whereas the alternative CBC with NCBP3 does not bind snRNA and associates only with mRNA thereby playing a role only in mRNA export. NCBP1/CBP80 is required for cell growth and viability (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RR205828