

Product datasheet for **RR205887**

Rab11fip4 (NM_001107023) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rab11fip4 (NM_001107023) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rab11fip4
Synonyms:	RGD1305063
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR205887 representing NM_001107023
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGCGGGCGGCCGGGCTGGGCGGGCGCCCCCGGGCTCTGCTGCGCTCCGTGCGCCGCTGCGCGAAG
TGTTTCGAGGTGTGCGGCCGCGACCCCGACGGCTTTCTACGCGTGGAGCGCGTCGCGGCCGCTCGGCCTGCG
CTTCGGCAAGGCGAGGAGGTGAAAAAATTGTGAAGTGTGGATCCCAATGACCTGGGAGAATCAAC
TTCAAGGATTTTTGCCGTGGAGTATTCGCCATGAAAGGGTGCAGGAGCTGCTGAAGGACGTGCTGTCGG
TGGAGAGTGCAGGGACCCTGCCATGTTCTCCAGACATCCAGACTGTGTGGAGCAGGGCAGCGACTTCTC
TGGCTCTACCTATGCTGATGGCAGCTCCTCCCAGGGAACCCGGCTTCTCCAGGAAGACGAGGAGGAA
GCGATGACATTGACCCTACCTGAGGGACCCAGGAGCTAGACATGGACAGCCCCATGGAGAGCAGCCAGG
GCCCGGAGGGCTCTGTCAGGAGTCTGTGGAGGAGAAGGAGCCTGAGCTGGGGGGCTGTTCTGCCAGA
GGACAAGTCCCTCGTCTCACTCCAGCCATGACCACTTCAGACCTTCCACTCACTCCACCGCCTCTCTC
ATCAGCAACGAGGAACAATTTGAAGACTATGGAGAGGGGACGAGCTGGACTGTGCTCCCAGCAGCCCTT
GCCCTGACGACGAGACCAGAACCAATGTCTACTCAGATCTAGGGTCATCTGTGTCTTCCAGTGCGGGACA
GACTCCCAGGAAGATGCGACATGCTTACAACAGTGACCTGCTAGATGTGTACTGCTCTCAGTGTGCAAG
AAGATCAACCTGCTGAATGACCTGGAGGCCAGGCTGAAAAACCTGAAGGCCAACAGTCCCAACAGGAAGA
TCTCAAGCACAGCTTTTGGACGGCAGCTCATGCACAGCAGTAACCTCAGCAGCAGCAATGGCAGCACCGA
GGACCTGTTCCGAGACAGCATCGACTCCTGTGACAATGACATCACAGAGAAGGTGAGCTTCTGGAAAAG
AAGGTAACGGAGCTGGAGAACGACAGCCTGGCGAGCGGGGGCTGAAGAGCAAGCTGAAGCTGGAGAACG
TGCAGCTGGTGCACAGGGTGCATGAGCTGGAGGAGATGGTGAAGGACCAGGAGACCAGCGGAGCAGGC
GCTGGAGGAGGAGGCCCGCGACACCGCGAGGTCTACTGCAAGCTGGAGCGGGAGAAGAGCACAGAGCTG
GAGCTGCTCAACACCAGGGTGCAGCAGTTAGAAGAAGAAAACACAGATCTGAGAACGACAGTGACACGAC
TCAAGTACAGACGGAGAAAAGTGGACGAGGAGCGGCAGCGCATGTCTGACCGGCTGGAGGACACCAGCCT
GCGGCTGAAGGACGAGATGGACCTTTACAAGCGCATGATGGACAAGCTGCGACAGAACCCTCGAGTTC
CAGAAGGAGCGGGAGGCGACGCAGGAGCTCATCGAAGACTTGGGAGAGAGCTGGAGCACCTGCAGATGT
ACAAGCTGGATTGCGAGCGGCCAGGCCGGGGCCGAGCTCATCCGGCCTTGGCGAGTTCAATGCCAGAGC
CCGTGAGGTGAACTTGAGCATGAGGTCAAGCGTCTAAGCAGGAGAATCACAACTGCGGGATCAGAAT
GATGACTTGAACGGGCAGATCCTGAGCCTCAGCCTCTATGAGGCAAAGAACCTTTTCCACCCAGACCA
AAGCCCAGTCTCTGGCTGCAGAAATCGACACAGCATCTCGAGATGAGCTAATGGAGGCCCTAAAGGAGCA
GGAGGAGATCAACTCCGACTGAGACAGTACATGGACAAGATTATCTGGCTATACTGGACCACAACCC
TCCATCTGGAGATCAAACAC

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR205887 representing NM_001107023
 Red=Cloning site Green=Tags(s)


MAGGAGWAGAPAAALLRSVRRLREVFEVCGRDPDGFLRVERVAALGLRFGQGEEVEKLVKCLDPNDLGRIN
 FKDFCRGVFAMKGCEELLKDVL SVESAGTLPSPDIPDCVEQGSDFSGSTYADGELLPREPGFFQEDEEE
 AMTLTLPEGPQELDMDSPMESSQGPEGSVRSPVEEKEPELGGLFLPEDKSLVLTAMTTSDLSTHSTASL
 ISNEEQFEDYGEADDVDCAPSSPCPDDETRNTVYSDLGSSVSSSAGQTPRKMRYHAYNSDLLDVYCSQCCK
 KINLLNDLEARKLNKANKSPNRKISSTAFGRQLMHSSNFSSSNGSTEDLFRDSIDSCDNDITEKVSFLEK
 KYTELENDSLASGGLKSKLKENVQLVHRVHELEEMVKDQETTAEQALEEEARRHREYVYCKLEREKSTEL
 ELLNTRVQQLEEEENTDLRTTVTRLKSQTEKLEERQRMDSRLEDTSRLKDEMDLYKRMMDKLRQNRLEF
 QKEREATQELIEDLRRELEHLQMYKLD CERPGRRSSSGLGEFNARAREVELEHEVKRLKQENHKL RDQN
 DDLNGQILSLSLYEAKNLFATQTKAQS LA AEIDTASRDELMEALKEQEEINFRLRQYMDKIILAILDHNP
 SILEIKH


SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:


ACCN: NM_001107023

ORF Size: 1911 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001107023.1](#), [NP_001100493.1](#)

RefSeq Size: 3127 bp
RefSeq ORF: 1914 bp
Locus ID: 303337
Cytogenetics: 10q25
MW: 72.1 kDa