

Product datasheet for RR205898L4

Psmc5 (NM_031149) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Psmc5 (NM_031149) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Psmc5
Synonyms:	Sug1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205898).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 81827
UniProt ID: [P62198](#)
Cytogenetics: 10q32.1
Gene Summary: may play a role in ATP-dependent RNA/DNA unwinding [RGD, Feb 2006]