

Product datasheet for **RR205962L4V**

Sox10 (NM_019193) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sox10 (NM_019193) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sox10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019193
ORF Size:	1398 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205962).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019193.2 , NP_062066.1
RefSeq Size:	3029 bp
RefSeq ORF:	1401 bp
Locus ID:	29361
UniProt ID:	O55170
Cytogenetics:	7q34
Gene Summary:	transcription factor associated with the sex determining factor SRY; may be important for glial cell development and cell specificity [RGD, Feb 2006]



[View online »](#)