

Product datasheet for **RR205968L3V**

Rtn4rl2 (NM_181380) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rtn4rl2 (NM_181380) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rtn4rl2
Synonyms:	Ngrh1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_181380
ORF Size:	1260 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205968).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181380.2 , NP_852045.1
RefSeq Size:	1304 bp
RefSeq ORF:	1263 bp
Locus ID:	311169
UniProt ID:	Q80WD1
Cytogenetics:	3q24
Gene Summary:	member of the reticulon 4 receptor family of proteins that may be involved in regulation of axonal regeneration [RGD, Feb 2006]



[View online »](#)