

## Product datasheet for **RR205990L3V**

### Ube2z (NM\_001037643) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2z (NM_001037643) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2z
Synonyms:	RGD1308347
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001037643
ORF Size:	1068 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR205990).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001037643.2</a> , <a href="#">NP_001032732.2</a>
RefSeq Size:	2874 bp
RefSeq ORF:	1071 bp
Locus ID:	303478
UniProt ID:	<a href="#">Q3B7D1</a>
Cytogenetics:	10q26
Gene Summary:	Catalyzes the covalent attachment of ubiquitin to other proteins. Specific substrate for UBA6, not charged with ubiquitin by UBE1. May be involved in apoptosis regulation. [UniProtKB/Swiss-Prot Function]



[View online »](#)