

Product datasheet for **RR205994**

Tmem259 (NM_001100886) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem259 (NM_001100886) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tmem259
Synonyms:	Mbrl; RGD1311136
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR205994 representing NM_001100886
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCGGAGCACGCGCAGCCCCAGGGCCCGGGCCAACGGCGCGGGGTGGCGCGCGGCCCGTGC
 GCGGCCGCGCGGCCCAATCTCAACCCCAACCGCTCATCAAGTGCGCGACCGCTCTCCACGCGCT
 CTTCTTCAAGATGGCCGTCACCTACTCACGCTCTTCCACCTGCTTCCGCGCCTCTTTGAGTCTTC
 GTTCTGCTCAAGGCCCTGTTGTGCTTTTCGCTCTCGCTACATACACATCGTCTTCTCCGGTCCCCA
 TCAACTGCTTGGAGCATGTTGAGACCGATGGCCACGGGAGGGTGTCTGCGGGTGGAGGTGCGCCAA
 CTCCAGCCGGCACCAGGTGATCTGCAGTTCTGTGATGGGGCCCTCGGTGGCCTAGACCTGGAGCCCGG
 GGCTGGAGCTGGAGGAGGAGGAGCTTCCGTGGAGATGTTACCAACAGCTCCATCAAGTTGAGCTGG
 ACATTGAGCCCAAGGTGTTCAAGCCACCAGGGGCGCCGATGCCCTGAATGACAGCCAGGACTTCCCTT
 CCCAGAGACACCAGCCAAAGCGTGGCCACAGGATGAGTACATTGTGGAGTACTCGCTGGAATATGGCTC
 CTGCGACTGTGCAAGCCACGCGCCAGCGCCTGAGCATCCCTGTATGTTAGTACCCCTAGACCCACGC
 GGGAGCAGTGCTTCGGGACCGCTTACGCCCTGCTGCTGGACGAGTTCCTGGGCTACGACGACGCTCT
 CATGTCTAGTGTGAAGGGCTGGCAGAGAATGAGGAGAACAAGGCTTCTTGGAGAAATGGTCTCTGGG
 GAGCACTACCGCTTTGTGAGCATGTGGATGGCAGCAGCATCTACCTGGCAGCCTTTGTATCATGGTCA
 TCTTTACCCTCAGTGTGCTCATGCTGTTGCGATACTCGCACCACCAGATCTTCGCTTTCATCGTGGACCT
 GCTGCAGATGCTGGAGATGAACATGGCCATCGCCTTCCCGCAGCGCCCTTGTGACCGTATCCTGGCT
 CTGCTCGGGATGGAGCCATCATGTCTGAGTCTTCAACGACACCACCAGCCCTTCTACATCATCTCA
 TCGTGTGGCTGGCCGACCAGTACGATGCCATCTGCTGCCACCAACACCAGCAAGCAGACTGAGCTGAG
 GTTCTTCTACCTTACCACTTCGCTTCTACGCTACCACTACCGCTTTAACGGCCAGTACAGCAGCTTG
 GCCCTGGTACGTCCTGGCTTTCATCCAGCATTCCATGATCTACTTCTTCCACCACTATGAGTTGCCCG
 CCATCTGCGAGCAGATCCGAATCCAGGAGATGCTGTTGCAGACACCCCGCTGGGCCCTGGGACCCAC
 GGCGCTGCCTGACGACCTCAACAATAACTCTGGCTCCCCTGCCACCCAGATCCCAGTCTCCCCTCGCG
 CTGGGCCCCAGCTCCAGCCCCGCGCCACTGGCGGAGCATCTGGGCTGGCTCACTGGGCGCTGGGGCT
 CAGTATCTGGCAGTGACCTAGGTTGGGTGGCCGAGACCGCCATCATCTCGGACGCGTCTTCTCTGTC
 GGGGCTGAGCGCCTCACTGCTGGAGCGCGCCAGCAGCCCCAGTCCCCGGACAGCGCAGCAGCTGAC
 CCTGGGCTCTCTGGAGGACGACCTGCCCTGCCGGTCC

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTAA

Protein Sequence:

>RR205994 representing NM_001100886
 Red=Cloning site Green=Tags(s)


MSEHAAAPGPGPNGGGGGAAPVRGPRGNLNPPLINVRDLFHALFFKMAVTSRLLFPFAFRRLFEFF
 VLLKALFVLFVLAYIHIVFSRSPINCLEHVRDRWPREGVLRVEVRHNSSRAPVILQFCDGGLGGLDEPG
 GLELEEEELSVEMFTNSSIKFELDIEPKVFKPPGADALNDSQDFPFPETPAKAWPQDEYIVEYSLEYGF
 LRLSQATRQRLSIPVMVVTLDPTREQCFGDRFSRLLLLDEFLLGYDDVLMSSVKGLAENEENKGFLLRVVSG
 EHYRFVSMWMARTSYLAAFVIMVIFTLVSMMLLRYSHHQIFVFIVDLLQMLEMNAIAFPAAPLLTVILA
 LVGMEAIMSEFFNDTTAFYIILIVWLADQYDAICCHTNTSKRHWRFFLYYHFAYHYRFNGQYSSL
 ALVTSWLFIQHSMIYFFHHYELPAILQQIRIQEMLLQTPPLGPGTPTALPDDLNNNSGSPATPDPSPPLA
 LGPSSSPAPTGGASGPGSLGAGASVSGSDLGWAETAIIISDASFLSGLSASLLERRPAAPSAPDSARP
 PGASLEDAPAPAGS

SGPTRTRPLEQKLI**SEEDLAANDILDYKDDDDKV**


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001100886
 ORF Size: 1722 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001100886.1 , NP_001094356.1
RefSeq Size:	2100 bp
RefSeq ORF:	1725 bp
Locus ID:	299608
Cytogenetics:	7q11
MW:	63.3 kDa