

Product datasheet for **RR206010**

Ccl4 (NM_053858) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccl4 (NM_053858) Rat Tagged ORF Clone

Tag: Myc-DDK

Symbol: Ccl4

Synonyms: Mip1-b; Scya4

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Cell Selection: Neomycin

ORF Nucleotide Sequence: >RR206010 representing NM_053858
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAAGCTCTGCGTGTCTGCCTTCTCTCTCTCTCTCTCTGCTTGTGGCCGCTTCTGCGATTCACTGCTGTCAG
CACCAATAGGCTCTGACCCTCCCACTTCCTGCTGCTTCTTTACACCTCCCGAAGATTCATCGAACTT
TGTGATGGATTACTATGAGACCAGCAGCCTTTGCTCCAGCCAGCTGTGGTATTTCTGACCAAAAAGGC
AGGCAAATCTGTGCCAGCCAGCGAGCCCTGGGTCAATGAGTATGTGAATGACTTGGAGTTGAAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR206010 representing NM_053858
Red=Cloning site Green=Tags(s)

MKLCVSAFSLLLLVAAFCDVLSAPIGSDPPTSCCFSYTSRKIHRNFVMDYYETSSLCSQPAVVFLTKKG
RQICADPSEPWNVEYVNDLELN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_053858

ORF Size: 276 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_053858.1</u> , <u>NP_446310.1</u>
RefSeq Size:	279 bp
RefSeq ORF:	279 bp
Locus ID:	116637
UniProt ID:	<u>P50230</u>
Cytogenetics:	10q26
MW:	10.2 kDa
Gene Summary:	a cytokine that is up-regulated during the inflammatory response; may be involved in intrapulmonary recruitment of neutrophils [RGD, Feb 2006]