

Product datasheet for **RR206016**

Fgg (NM_012559) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fgg (NM_012559) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fgg
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR206016 representing NM_012559 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATTGGTCCTTGCAACTCCGGAGTTTCATTCTATGCTGGCGCTTTACTACTTTTCCAACGGGCC
TGGCACAGTACACTGCTACTAGAGATAACTGTTGCATCCTAGATGAGAGATTCGGTAGTTACTGCCAAC
CACTTGTGGCATCTCAGACTTCTGAATTCCTACCAACCGACGTGGACACTGACCTCCAGACTCTGGAA
AACATCTTACAACGAGCTGAAAACAGAACCACAGAAGCCAAGGAACTGATTAAGCCATCCAGGTTTACT
ACAACCCAGACCAACCCCAAGCCGGGTATGATAGAGGGTGCCACTCAGAAGTCTAAGAAGATGGTAGA
AGAAATCTTGAAGTACGAAGCACTGTTACTAACCATGAGTCGAGTATCCGGATTTACAGGACATCTAT
ACTTCAAATAAGCAGAAGATCACGAACCTGAAACAGAAGGTGGCCAGCTTGAAGCTCAGTGCCAGGAGC
CTTGCAAGGACTCTGTGCGAATCCATGACACAACCGAAAAGATTGTCAGGATATCGCCAACAAGGGTGC
CAAAGAAAGTGGACTTTACTTCATTTCGGCCTTTGAAAGCCACGCAGCAGTTCTTAGTGTACTGTGAAATC
GACGGATCTGAAATGGCTGGACCGTATTGCAGAAGAGGCTTGATGGCAGTGTGGATTTCAAGAAGAACT
GGATTCAGTATAAAGAAGGATTTGGACACCTGTCACCTACTGGCACCACAGAGTTTTGGCTGGGCAATGA
GAAGATTCATCTGATAAGCATGCAGTCCACCATCCCATACGCACTGAGAATACAGCTCAAAGACTGGAGT
GGCAGGACCAGCACCGCGGACTATGCCATGTTCCAGGGTGGTCCCTGAATCCGACAAATACCGCTGACCT
ATGCCTACTTCATTGGCGGAGATGCCGGGGATGCCTTCGACGGCTACGATTTTGGTGATGATCCCAGTGA
CAAGTTTTTACATCCCACAACGGCATGCACCTTCACTGACCTGGGACAATGACAACGACAAGTTTCAAGGC
AACTGTGCTGAGCAGGATGGATCTGGCTGGTGGATGAACAAGTGTACGCTGGCCACCTCAATGGAGTTT
ATTACCAAGGTGGCACTTACTCCAAGTCATCTACTCCTAACGGTTATGACAATGGCATTATTTGGGCCAC
CTGGAAAACCCGCTGGTATTCCATGAAGGAAACCACATGAAGATAATTCCTTCAACAGACTCTCCATT
GGAGATGGCAGCAACATCACATGGGAGGATCCAAACAGGTTGGAGACATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR206016 representing NM_012559
 Red=Cloning site Green=Tags(s)


MNWSLQLRSFILCWALLLLSPTGLAQYTATRDNCCILDERFGSYCPTTCGISDFLNSYQTDVDTDLQTL
 NILQRAENRTTEAKELIKAIQVYYPDQPPKPGMIEGATQKSKMVVEILKYEALLLTHESSIRYLQDIY
 TSNKQKITNLKQKVAQLEAQCQEPCKDSVRIHDTTGKDCQDIANKGAKESGLYFIRPLKATQQFLVYCEI
 DSGSGNGWTVLQKRLDGSVDFKKNWIQYKEGFHLSPTGTTEFWLGNEKIHLISMQSTIPYALRIQLKDW
 SGRSTADYAMFRVGPESDKYRLTYAYFIGGDAGDAFDGYDFGDDPSDKFFTSHNGMHFSTWDNDNDKFE
 G NCAEQDGSWWMNKCHAGHLNGVYYQGGTYSKSSTPNGYDNGI IWATWKTRWYSMKETTMKIIPNRLSI
 GDGQQHHMGGSKQVQDM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_012559
ORF Size:	1311 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012559.2 , NP_036691.2
RefSeq Size:	1554 bp
RefSeq ORF:	1314 bp
Locus ID:	24367
UniProt ID:	P02680
Cytogenetics:	2q34
MW:	49.7 kDa
Gene Summary:	combines with alpha and beta polypeptides to form the mature fibrinogen molecule; plays a role in inflammatory response [RGD, Feb 2006]