

Product datasheet for **RR206049L4V**

B4galt4 (NM_001012018) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	B4galt4 (NM_001012018) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	B4galt4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001012018
ORF Size:	1032 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206049).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001012018.1 , NP_001012018.1
RefSeq Size:	2106 bp
RefSeq ORF:	1035 bp
Locus ID:	303923
UniProt ID:	Q66HH1
Cytogenetics:	11q21
Gene Summary:	Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.[UniProtKB/Swiss-Prot Function]



[View online »](#)