

## Product datasheet for **RR206060L3V**

### Hsp90ab1 (NM\_001004082) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hsp90ab1 (NM_001004082) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hsp90ab1
Synonyms:	HSP90-BETA; HSP90B; HSPC2; Hspcb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001004082
ORF Size:	2172 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206060).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001004082.3</a> , <a href="#">NP_001004082.3</a>
RefSeq Size:	2338 bp
RefSeq ORF:	2175 bp
Locus ID:	301252
UniProt ID:	<a href="#">P34058</a>
Cytogenetics:	9q12
Gene Summary:	heat shock protein; displays chaperone activity; binds to peptides and facilitates their appropriate transport and/or secondary structure [RGD, Feb 2006]



[View online »](#)