

Product datasheet for **RR206061**

Cpa4 (NM_001109346) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cpa4 (NM_001109346) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cpa4
Synonyms: RGD1563619
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR206061 representing NM_001109346
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGTGGTTACTGTTCTTTGGGGCCCTTATCGGGACAGGCATCTGCGGCCGAGACAAGTTCTTTGGGG
ACCAAGTTTTTAGGATTAATGTCAGAAATGGAGACGAGATCAGGAACTGACCGAGCTAGTGAATTCGGA
CCACTTTAAGCTCAATGTCTGAAATCTCCCTCCACTTTTGATCGGCCTGTGGATATTCTCGTCCCTTCT
GTCAGCTTGCTGCCAGTCAAGTCTTCCCTGAAAGTCCCAGGGTTTAGACTACTCAGTGACAATTGACAACC
TACAGGCTCTCTTAGATATTGAAGATGAAGAAATGCAACACAATGAAGGGAGAGAACGGAGTGGTGA
CAACTACGGATCTTACCATTCCCGGAGGCTATTTACCACGAGATGGACAGTATTGCCACAGACTTCCCT
GACCTGGCAAGCCGAGTGAAGATCGGAGAGACATTCGAAAAGCGGCCATGTATGTTCTCAAGTTCAGCA
CGGGAGGAGGAGGCAAGAAACGGCCAGCCATTTGGCTGAATGCAGGCATCCATGCCCGTGAAGTGGATCTC
ACAGGCCACTGCCATCTGGACAGCAAGGAAGATTGTGACTGATTATCAGAAGGACCCAGCTGTCACTCC
ATTTTAGAGAAAATGGATATTTTCTATTGCCTGTGGCCAATCCTGACGGATATGTGTACACAAAAAC
AAAACCGATTGTGGAGGAAGACACGGTCCCGGAACCCAGGAAGCCGCTGTATTGGAGCTGATCCAACAC
AAACTGGAATGCTAGTTTTGCAGGAGAGGGGCCAGCAATAACCCCTGCTGAAGTATACCATGGCTCC
CACCCCAATTCTGAAGTGGAGGTGAAATCGGTGGTGGACTTTATTCAAAAGCATGGGAACCTTCAAATGCT
TCATTGACCTGCATAGCTACTCACAGCTGCTGATGATCCCTATGGTACTCGGTCAAGAAGGCCCCAGA
TGCCGAGGAGCTGGACGATGTGGCAGGAGTGCAGCCAAAGCCCTGGCTTCTCTCTCGGGAACCAAGTAC
CAAGTGGGACCGACTTGTACCACCGTCTATCCAGCTAGCGGGAGCAGCGTTGACTGGGCATACGACAACG
GCATCAAGTATGCTTTCACATTTGAGCTGAGAGACTGGGCACTATGGCTTCTCTTGGCAGCCAGCCA
GATCATCCCCTGCGGAGGAGACTGGCAAGGGCTGAAGGTCATCATGGAACACGTTAGGGACAACCTC
TAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR206061 representing NM_001109346
 Red=Cloning site Green=Tags(s)


MKWLFFGALIGTGICGRDKFFGDQVFRINVRNGDEIRKLTENVSDHFKLNVWKPSTFDRPVDILVPS
 VSLLPVKSFLLKSQLDYSVTIDNLQALLDIEDEEMQHNEGRERSGDFNYGSYHSPEAIYHEMDSIATDFP
 DLASRVKIGETFEKRPYVLFKSTGGGKKRPAIWLNAGIHAREWISQATAIWTARKIVTDYQKDPAVTS
 ILEKMDIFLLPVANPDGYVYTQNQNRLLWRKTRSRNPGSRCIGADPNRNWNASFAGEGASNNPCSEVYHGS
 HPNSEVEVKSVDVFIQKHGNFKCFIDLHSYSQLLMYPYGYSVKKAPDAEELDDVARSAAKALASLSGTKY
 QVGPTCTTVYPASGSSVDWAYDNGIKYAFTELRDTHYGFLLPASQIIPTAEETWQGLKVI MEHVRDNL
 Y

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_001109346

ORF Size:	1263 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001109346.2 , NP_001102816.1
RefSeq Size:	2781 bp
RefSeq ORF:	1266 bp
Locus ID:	502736
Cytogenetics:	4q22
MW:	47.4 kDa