

Product datasheet for **RR206075L4V**

Plac1 (NM_001024894) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Plac1 (NM_001024894) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Plac1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001024894
ORF Size:	519 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206075).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024894.1 , NP_001020065.1
RefSeq Size:	1103 bp
RefSeq ORF:	522 bp
Locus ID:	317316
UniProt ID:	Q4V7E2
Cytogenetics:	Xq37
Gene Summary:	May play a role in placental development.[UniProtKB/Swiss-Prot Function]



[View online »](#)