

Product datasheet for **RR206090L4V**

Slc6a12 (NM_017335) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc6a12 (NM_017335) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc6a12
Synonyms:	BGT1; Gat1; RNU28927; VGAT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017335
ORF Size:	1884 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206090).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017335.1 , NP_059031.1
RefSeq Size:	2561 bp
RefSeq ORF:	1887 bp
Locus ID:	50676
UniProt ID:	P48056
Cytogenetics:	4q42
Gene Summary:	sodium/chloride-dependent neurotransmitter transport protein; expressed in exocytotic vesicles in GABAergic cells; associated with exocysts [RGD, Feb 2006]



[View online »](#)