

Product datasheet for **RR206120**

Suclg2 (NM_001100750) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Suclg2 (NM_001100750) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Suclg2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR206120 representing NM_001100750
Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC**GGCGC**
GCCC

ATGGCGTCCGCGGTGCCATCGCGGCCAGGCCGGGAAGCTTCTGCGAGAGCGATCGCTGCGGCCCCGAC
TCTTGGCGGTCAAGTCCCAGGCAGTCACTTAACCCCAAGGAGATGGCTGAACCTGCAGGAATACCAGAG
CAAGAAGTCATGTCAGAGCACGGAGTGAGAGTGCAGAGGTTCTTTGTCGCCAGCACTGCTAAGGAGGCT
CAGGAAGCAGCCAAGCGACTGAATGCAAAGGAAATTGTTTTAAAAGCCAGATCTTAGCTGGAGGAAGAG
GAAAAGGTGTCTTCGATAGCGGTTTAAAAGGAGGTGTTCAATTAACAAAAGACCCATAAAGTGGTGGGACA
GCTGGCTCAGCAGATGATTGGTTATAATCTAGCAACAAAACAACTCCAAAAGAAGGTGTGAAAAGTTAAC
AAGGTGATGGTTGCTGAAGCCCTGGATATTTCTAGAGAAACGTAAGTGGCCATTCTGATGGACCGGTCCC
ATAATGGCCCTGTGTTAGTAGGCAGTCCCAGGGGGCGTCGATATTGAAGAAGTGGCAGCTTCCAGTCC
AGAAGTCACTTTAAGGAGCAGATTGACATCTTTGAAGGATAAAGGATAGCCAGGCCAGCGCATGGCG
GAAAATCTGGGCTTCTCGGGTCTTTGAAAACCAGGCTGCAGATCAGATTAAGGCTGTACCATCTCT
TCCTGAAGATTGACGCCACTCAGGTAGAAGTGAACCCCTTCGGTAAAACCCAGAAAGGCAAGTTGTCTG
TTTTGATGCCAAGATAAACTTTGATGACAATGCTGAATTCGACAAAAGGACATATTTGCTATGGACGAC
AAGTCAGAGAATGAACCCATTGAAAATGAAGCTGCCAGATATGATCTCAAATATATAGGTCTAGACGGGA
ACATCGCCTGCTTTGTGAATGGTGTGGATTGGCCATGGCTACTTGTGACATCATTTTCTGAAATGGAGG
GAAGCCAGCTAAGTCTTGGACCTCGGAGGTGGTGTAAAGGAGTCCCAAGTCTATCAAGCATTCAAAGT
CTCACATCTGATCCCAAGGTGGAAGCCATTCTTGTCAATATCTTTGGTGGGATCGTCAACTGTGCCATCA
TTGCCAACGGAATCACCAAAGCCTGCCAGGAGCTAGAGCTCAAGGTGCCACTGGTAGTCCGTCTTGAAGG
AACTAATGTCCAGGAGGCTCAGAACATCTCAAAGCAGCGGACTCCCCATTACGTACAGTGTGGACCTG
GAGGATGCAGCCAAGAAAGCTGTTGCCAGTGTGGCCAAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR206120 representing NM_001100750
 Red=Cloning site Green=Tags(s)


MASAVAIAAQAGKLLRERSLRPRLAVRSQAGHLTPRRWLNQEQYQSKLMSEHGVRVQRFFVASTAKEA
 QEAAKRLNAKEIVLKAQILAGGRGKGVFDSGLKGGVHLTKDPKVVGQLAQQMIGYNLATKQTPKEGVKVN
 KVMVAEALDISRETYLAILMDRSHNGPVLVGGSPQGGVDIEEVAASPELIFKEQIDIFEGIKDSQAQRMA
 ENLGFLGSLKNQAADQIKKLYHLFLKIDATQVEVNPFGETPEGQVVCDAKINFDDNAEFRQKIDIFAMDD
 KSENEPIEENEAARYDLKYIGLDGNIACFVNGAGLAMATCDIIFLNGGKPANFLDLGGGVKESQVYQAFKL
 LTSDPKVEAILVNIFGGIVNCAIIANGITKACQELELKVPLVVRLEGTNVQEAQNILKSSGLPITSAVDL
 EDAAKKAVASVAKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

AscI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001100750
ORF Size:	1302 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001100750.1 , NP_001094220.1
RefSeq Size:	2263 bp
RefSeq ORF:	1305 bp
Locus ID:	362404
Cytogenetics:	4q34
MW:	47 kDa