

Product datasheet for **RR206171L3V**

Apba1 (NM_031779) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Apba1 (NM_031779) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Apba1
Synonyms:	Mint1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031779
ORF Size:	2517 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206171).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031779.2 , NP_113967.1
RefSeq Size:	3681 bp
RefSeq ORF:	2520 bp
Locus ID:	83589
UniProt ID:	O35430
Cytogenetics:	1q51
Gene Summary:	component of multimeric complex containing Munc18-1 and syntaxin; may mediate synaptic vesicle docking/fusion [RGD, Feb 2006]



[View online »](#)