

Product datasheet for **RR206188L3V**

Rab31 (NM_145094) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rab31 (NM_145094) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rab31
Synonyms:	Rab0
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145094
ORF Size:	585 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206188).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145094.2 , NP_659562.2
RefSeq Size:	3335 bp
RefSeq ORF:	588 bp
Locus ID:	246324
UniProt ID:	Q6GQP4
Cytogenetics:	9q37
Gene Summary:	an oligodendrocyte Rab protein [RGD, Feb 2006]



[View online »](#)