

Product datasheet for **RR206221L4V**

Igsf11 (NM_001013120) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Igsf11 (NM_001013120) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Igsf11
Synonyms:	MGC94879
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001013120
ORF Size:	1284 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206221).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001013120.1 , NP_001013138.1
RefSeq Size:	3675 bp
RefSeq ORF:	1287 bp
Locus ID:	303926
UniProt ID:	Q5U2P2
Cytogenetics:	11q21
Gene Summary:	Functions as a cell adhesion molecule through homophilic interaction. Stimulates cell growth (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)