

Product datasheet for **RR206225**

Sephs1 (NM_001104630) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sephs1 (NM_001104630) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sephs1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR206225 representing NM_001104630 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTACTCGGGAGTCCTTTAACCCGGAACTTATGAATTGGACAAGAGCTTCCGGCTAACCGGTTTA
CTGAACTGAAAGGCACAGGCTGCAAAGTGCCCAAGATGTCTTGCAAGAAATTGCTGGAATCTTTACAAGA
GAACCACTTCCAAGAAGATGAGCAGTTTCTGGGAGCGTTATGCCACGACTTGGCATTGGGATGGATACT
TGCCTATTCTCTGAGGCACGGTGGTCTCTCCTTGGTTCAAACCACAGATTACATTTATCTATCGTGG
ACGACCCCTTACATGATGGCAGGATAGCGTGTGCCAATGTCCTCAGTGACCTCTATGCGATGGGTGCAC
AGAATGTGACAATATGCTAATGCTCCTTGAGTCAAGTAATAAAATGACTGACCGGAAAGGATAAAGTG
ATGCCGCTAATTACAGGGTTTTAAAGATGCGGCCGAGGAAGCGGGACCTCTGTAACGGGGGCCAAA
CGGTATTGAACCCCTGGATTGTTCTGGGAGGAGTCGCCACAAGTCTGCCAGCCCAATGAATTTATCAT
GCCAGACAACGCGGTGCCGGGGACGTGCTGGTGTGACCAAACCCCTGGGGACGCAGGTTGCAGTGGCT
GTGCACCAGTGGCTGGATATTCTGAAAAATGGAATAAAATTAAGCTCGTGGTCACCCAAGAAGACGTAG
AGCTGGCATAACCAAGAGGCAATGGTGAACATGGCGGGCTCAACAGGACAGCTGCAGGCCTCATGCACAC
GTTCAATGCTCATGCAGCCACTGATATCACAGGCTTCGGGATTCTGGCCATGCACAGAACCTGGCCAAG
CAACAGAGGAACGAGGTGTCGTTTGTGATTCACAATCTTCCCGTGTAGCCAAGATGGCCGCTGTGAGCA
AAGCCTGCGGAAACATGTTTCGGCCTAATGCATGGGACCTGCCAGAGACGTGAGGAGGCCTGCTAATCTG
CCTACCGCGTGAGCAAGCAGCTCGGTTCTGTGACAGAGATAAAGTCCCCAAATATGGCGAAGGCCACCA
GCATGGATTATTGGGATTGTAGAGAAGGGCAACCCGACAGCCAGAATCATCGACAAACCTCGGATTATTG
AAGTCGCACCTCAAGTAGCCACGCAAAACGTGAATCCCACACCTGGTGCCACCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR206225 representing NM_001104630
 Red=Cloning site Green=Tags(s)

MSTRESFNPETYELDKSFRLTRFTELKGTGCKVPQDVLQKLLSLQENHFQEDEQFLGAVMPRLGIGMDT
 CVIPLRHGGLSLVQTTDYIYPIVDDPYMMGRIACANVLSLYAMGVTECDNMLMLLGVSNKMTDRERDKV
 MPLIIQGFKDAEEAGTSVTGGQTVLNPWIVLGGVATTVCQPNFIMPDNAVPGDVLVLTPLGTQVAVA
 VHQWLDIPEKWINKIKLVVTQEDVELAYQEAMVNMARLNRTAAGLMHTFNAHAATDITGFGLGHAQNLAK
 QQRNEVSFVIHNLPLAKMAAVSKACGNMFLMHGTCPE TSGGLLIICLPREQAARFCAEIKSPKYGEGHQ
 AWIIGIVEKGNRTARIIDKPRIIEVAPQVATQNVNPTPGATS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001104630

ORF Size: 1176 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001104630.1](#), [NP_001098100.1](#)

RefSeq Size: 5263 bp


RefSeq ORF: 1179 bp

Locus ID: 291314

Cytogenetics: 17q12.3

MW: 42.9 kDa

Product images:



Circular map for RR206225