

Product datasheet for **RR206232**

Zfp593 (NM_001106689) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Zfp593 (NM_001106689) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Zfp593
Synonyms: Znf593
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR206232 representing NM_001106689
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGGGTCGCTCCCGCGGACCGCGCGCACCGAGCACATTCCTTAGCCCGCCAGATGAAGGCCAAGAAGC
GGCGGCCGACCTGGATGAGATCCATCGCGAGCTGAGGCCGAGGCTCTCCGCGTCCAAAGCGGGAGCG
GGACGCAGAGCCGGATCCGGACCTGCCAGGGGGCGGCTGCATCGCTGTCTGGCCTGCGCGAGGTACTTC
ATCGATTCTGCCAACCTGAAGACCACTCCGATCCAAAGACCACAAGAAAAGGCTGAAGCAGCTGAGTG
TCGAACCTACAGTCAGGAAGAGGCTGAGAGGGCCGCGGCATGGGCTCATATGTGCAACCTCAGCGACT
GGCGTGCCACAGAAGTGTCCACTGAGATCCCTGAGATGGACACGTCCACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR206232 representing NM_001106689
Red=Cloning site Green=Tags(s)

MGRSRRTGAHRAHSLARQMKAKKRRPDLDEIHRELRPQALPRPKRERDAEPDPLPGGGLHRCLACARYF
 IDSANLKTHFRSKDHKKRLKQLSVEPYSQEEAERAAGMGSYVQPQRLGVPTEVSTEIPEMDTST


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM 001106689

ORF Size: 402 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106689.1](#), [NP_001100159.1](#)


RefSeq Size: 2010 bp

RefSeq ORF: 405 bp

Locus ID: 298546

Cytogenetics: 5q36
MW: 15.3 kDa

Product images:



Circular map for RR206232