

Product datasheet for **RR206233L3V**

Map4 (NM_001024278) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Map4 (NM_001024278) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Map4
Synonyms:	Mtap4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001024278
ORF Size:	3171 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206233).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024278.1 , NP_001019449.1
RefSeq Size:	5535 bp
RefSeq ORF:	3174 bp
Locus ID:	367171
UniProt ID:	Q5M7W5
Cytogenetics:	8q32
Gene Summary:	Non-neuronal microtubule-associated protein. Promotes microtubule assembly (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)