

## Product datasheet for **RR206246**

### Rpl10a (NM\_031065) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rpl10a (NM\_031065) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Rpl10a  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR206246 representing NM\_031065  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGAGCAGCAAAGTTTCACGCGACACCCTCTACGAGGCGGTGCGGGAAGTCTGCACGGGAACCAGCGTA  
AGCGCCGCAAGTTTCTGGAGACGGTGGAGCTGCAGATCAGCCTGAAGAACTACGACCCTCAGAAGGACAA  
ACGCTTCTCGGGCACCGTCAGGCTTAAGTCCACCCTCGCCCCAAGTTCTCCGTGCGTCTGGGGGAC  
CAGCAGCACTGTGATGAAGCCAAGGCCGTGGATATCCCCACATGGACATCGAGGCGCTCAAGAAGCTTA  
ACAAGAACAAGAAGTTGGTCAAAAAGCTGGCCAAGAAGTATGATGCCTTTTTGGCCTCCGAGTCTCTGAT  
CAAGCAAATCCCACGTATCCTGGGCCAGGCTAAACAAGGCGGCAAGTTCCCCTCCCTGCTGACACAC  
AATGAAAACATGGTGGCCAAAGTGGATGAGGTGAAATCCACAATCAAGTCCAGATGAAGAAGGTGCTGT  
GTCTGGCCGTGGCGTCGGCCACGTGAAGATGACCGATGACGAGCTAGTCTACAACATCCACCTGGCTGT  
CAACTTCTGGTGTCTTACTCAAGAAAAGTGGCAGAAGCTCCGGGCTCTGTACATCAAGAGCACCATG  
GGCAAGCCCCAGCGTCTGTAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR206246 representing NM\_031065  
Red=Cloning site Green=Tags(s)

MSSKVS RD TLYEAVREVLHGNQRKRRK FLETVELQISLKNYDPQDKRFS GTVRLKSTPRPKFSVCVLGD  
QQHCDEAKAVDIPHMDIEALKKLNKKNKLYKKLAKKYDAFLASESLIKQIPRILGPGLNKAGKFP SLLTH  
NENMVAKVDEVKSTIKFQMKKVLCLAVAVGHVKMTDDEL VYNIHLAVNFLV SLLKKNWQNV RALYIKSTM  
GKPQRLY


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_031065

ORF Size: 651 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_031065.1</a></u> , <u><a href="#">NP_112327.1</a></u>
<b>RefSeq Size:</b>	710 bp
<b>RefSeq ORF:</b>	654 bp
<b>Locus ID:</b>	81729
<b>UniProt ID:</b>	<u><a href="#">P62907</a></u>
<b>Cytogenetics:</b>	20p12
<b>MW:</b>	24.8 kDa
<b>Gene Summary:</b>	60S ribosomal subunit protein [RGD, Feb 2006]