

Product datasheet for **RR206260**

Fcho1 (NM_001106069) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fcho1 (NM_001106069) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fcho1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR206260 representing NM_001106069
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGGAAGGCTTCACTGTGCGCCCCGATACATCTCAGAACAAATATCCTTTGACTCTGCCGAGCCCCAC
 GCTTCTCCTCTAGTGACTCTGATTTTGATGATGAGGAGCCTAGAAAGTTCTATGTTGCGATCAAGCCTGC
 CCCC GCCCGACCATGGCCTGCAGCTCTGAAGCTGCAGCCGCCAGCTCAGGGCCACAGCTGGCAGCCTC
 ATCCTCCCGCCAGGCCAGGGGGCACCATGAAACGTCTCTCATCGCGGGACACTTCTGGGAAGCCACAGA
 GACCTCGGTGGCCCCACGCAGCAGCTGCACAGAGAAGCCGCTGGCCTCGGAGGAGCCACTGTCCAAGAG
 CCTCTTTGGCGCCGCTGGAGTCTGCCTTCGACCAGGAGACTTCACAGGCTCCAGCAGCCTGGGCTTC
 ACATCCAGCCCATCTCCTTTCTCATCTCATCTCCGAGAATGTGAAGACTCGGGCCTGGACTCTCCAT
 CTCTTGCTGCACCCGGTCCATCGCCTGAGTCTGGTCCCCGACCAGGCACCCGCAGAGCCCACCCAC
 CTGCAGGGCACAGCACCAGGAGGCCAGGGACCCAGTGCCTCGTGCCTCGCCAGGTCCCTGGGGCCT
 GAGGCAGGTGCAGGTGCAGACTCACTGACGCTGGCTGACCCACCAGGGAGGGCCTGGCTGCCACGCTGA
 GGAGACCACGGTCCAGGAAAGTGTCTGCCCGCTCACCAGGAGCAGTGGTGACCTGTGCCGTCGCTCAG
 CCCATCCCCACTGGGATCCTCTGCCCCACCATCCCCCCTGACCGGCCAGCTTCTCCGCCAGATGGGA
 CACGGTGTCTCCCGTGGCCCCAGCCAGTGGTCTGGGATCCCAGGACACCTTGCTGTGCCACAGCCT
 TCACTGAGTATGTCTATGCCTTTTTCCGTGGCCACAACCCAGTTGCCTGGCTCGAGTACCCGGGAGTT
 GACCATGACCTTCCCTGCGGGCATTGTACGTGTGTTCACTGGGACCCACCACCACTGTCTCAGCTTC
 CGGCTGGTGCATACGTCCCCATAGAACACTTACAGCCCAACGCTGACCTGATCTTCAGCGACCCCTCCC
 AGAGTGACCCGGAGACCAAGATTTCTGGCTGAACATGGCGGCATTGACGAGGCCCTGCAGCACCAGGC
 AGAGCAGAACCCACCGCCTCCTACTACAACCTGGTGTGTTGCGGTACCAAGTTCTCCCGCCTGGACCC
 AAGTCAATGCCCTGCAATGAGCGCCCACTGGCAGTGCGGGCCACACTCACGCGGGTTTCGGTGGAGT
 ATAGCTACCGTGCGGGTGCCACTGCTGTGTCCACGCCACTCACTAACGTCCAGATCCTGTGCCGTGGG
 GGAGCCGGTACCAGTGTGCGGTGCAGCTGCGGCCACCTGGAACACAGAGGAAAAGAGTTTACATGG
 AAGCTTCCAGATGTGCGTGAGGCAGGGGCTCGGGCCACCTGTCCGCCAGCTGGCAGCCGAGTCCGGGC
 CCAGCACCCCTAGCCCCGTGACCGCGCAATTCACCACTGAGGGCGCCACGCTATCCGGCCTGGACCTGGA
 GCTGCTGGGTGGCGGTACCGCATGTCCCTGGTGAAGCGGAGATTCACAGGGATGTACCTCACCAGCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR206260 representing NM_001106069
 Red=Cloning site Green=Tags(s)


MRKASLCAPIHLRTISFDSAEPFRFSSSDSDFDDEEPRKFYVRIKPAPARTMACSSEAAAAQLRATAGSL
 ILPPGPGGTMKRLSSRDTSKPKQRPSPRSPSSCTEKPLASEEPLSKSLFGPPLESFDFHEDFTGSSSLGF
 TSSPSPFSSSSPENVEDSGLDSPSLAAGPSPESWVPRPGTPQSPPTCRAQHQAARDPVPRVPSPGPWGP
 EAGAGADSLTLADPTREGLAATLRRPRSRKVSCLTRSSGDLCRSLSPSPLGSSAPTIPPDRPSFSAQMG
 HGVSRRGSPVVLGSQDTLPVATAFTEYVYAFFRGHNPSCLARVTGELTMTFPAIVRVFSGTPPPPVLSF
 RLVTSPVIEHLQPNADLIFSDPSQSDPETKDFWLNMAALTEALQHQAEQNPTASYYNLVLLRYQFSRPGP
 KSMPLQMSAHWQCGPTLTRVSVEYSYRAGATAVSTPLTNVQILLPVGEPVTSVRLQPAATWNTTEKRFTW
 KLPDVREAGSGHLSASWQPQSGPSTSPVTAQFTSEGATLSGLDLELLGGGYRMSLVKRRFTGMYLTSC

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001106069

ORF Size: 1680 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106069.1](#), [NP_001099539.1](#)

RefSeq Size: 3236 bp


RefSeq ORF: 1683 bp

Locus ID: 290639

Cytogenetics: 16p14

MW: 59.9 kDa

Product images:



Circular map for RR206260