

Product datasheet for **RR206272**

Dpt (NM_001105965) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dpt (NM_001105965) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dpt
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR206272 representing NM_001105965
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACCTCACTCTTCTGTGGGTGCTTCTGCCACTGGTGAGCACAGCCTGGGGCCAGTATGGTGACTACG
GGTACCATCCATACCAGCAGTACCAGGACTACAGTGACGACGGATGGGTGAACCTGAACCGTCAGGGCTT
CAGCTACCAGTGTCCCCACGGCAAGTGGTGGTAGCTGTGAGGAGCATCTTTAACAAAGAAGGAGGGCTT
GACAGACAGTGGAACTATGCCTGCATGCCCACTCCCCAAAGCCTCGGGGAGCCCACTGAGTGCTGGTGGG
AGGAGATCAACAGGGCTGGCATGGAATGGTACCAGAAAGTCTTAACAATGGATTGGTGGCGGGTTCCA
AAGCCGTTACTTCGAGTCAGTGTGGATCGAGAGTGGCAATTTTACTGCTGCCGGTATAGCAAGAGGTGC
CCATATTCCTGCTGGGATCGCCAGTGGAAAGTTCATAATGTGCCGGATGACTGACTATGACTGTGAATTCC
AAAACGTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR206272 representing NM_001105965
Red=Cloning site Green=Tags(s)

MDLTLLWVLLPLVSTAWGQYGDYGYHPYQQYQDYSDDGWVNLNRQGSYQCYPHGQVVAVRSIFNKKEGS
DRQWNYACMPTPQSLGEPTECWEEINRAGMEWYQKCSNGLVAGFQSRYPFESVLDREWQFYCCRYSKRC
PYSCWDRQWKFIMCRMTDYDCEFENV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001105965

ORF Size: 498 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105965.1, NP_001099435.1</u>
RefSeq Size:	1618 bp
RefSeq ORF:	501 bp
Locus ID:	289178
Cytogenetics:	13q23
MW:	20 kDa