

Product datasheet for **RR206274L3V**

Tnfrsf1b (NM_130426) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Tnfrsf1b
Synonyms:	Tnfr2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_130426
ORF Size:	1444 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RR206274).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_130426.4
RefSeq Size:	2924 bp
RefSeq ORF:	1425 bp
Locus ID:	156767
UniProt ID:	Q80WY6
Cytogenetics:	5q36



Gene Summary: mediates brain uptake of tumor necrosis factor (TNF)-alpha from the circulation [RGD, Feb 2006]