

Product datasheet for **RR206275**

Olr67 (NM_001001279) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olr67 (NM_001001279) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olr67
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR206275 representing NM_001001279 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTCTGTTTCAACAACCTCTGAGGTCAAGTATTTCTTCTAATTGGGATCCCAGGACTGGAATATGCC
ATGCATGGATCTCCATCCCCATCTGCCTCATGTACCTATTGCCATCGTGGGTAAGTGCATATCATTTT
TGTCATCAAGACAGAGCCTTCTCTTCATGAACCATGTACTATTCCTTACCATGTTGGCTGTCTCTGAT
ATGGGTTTGCTCTCTCTCCCTTCTACCATGCTTAAGATCTTCTGTTCAATGCCATGCAAATCTCAC
CCAATGCCTGTTTGTCAAGAGTTCTTTATCCATGTATTCAGTCCATGGAGTCTCTGTGCTTCTGAT
TATGTCTTTGGACGCTTCTTGCCATTCAACCCCTTAAGGTACAGTTCTATCCTCACTGGCAGTAGG
GTTGCTAAATAGGGTTAATTTTAGCCATCAGGAGCACAGCTTTGGTGCTTCCATTCCCTTTCACACTAA
GGAACTGAAATACTGTCAAAGAATCTGCTCTCTCACTTATACTGCCTTCATCAGGACGTGATGAACT
GGCCTGCTCTGATAACAAGATCAATTTCTATGTTTCTTTGCGGCTCTCTGTACTATGCTGGACTTT
GCCTTGATTCTTATGTCTTACATGCTGATCTTGAAGACAGTTCTCAGTATTGCATCTCTGGTAGAAAGGC
TAAAGCCCTAAATACCTGTGTCTCCACATCTGTGCTGTGCTCATCTTCTATGTCCCATCATCACCTT
GGCTGCCATACATCGCTTTGCCAAGCACAAAAGCCCCCTACTTGAATCCTTATTGCAGATATGTTTTG
TTGGTGCCACCCCTGATGAACCCATTGTGTACTGTATAAAGACTCGACAAATCCGGGAGAAGGTCTTGG
GAAAGTTGGTAACTTATGTGTGAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RR206275 representing NM_001001279

Red=Cloning site Green=Tags(s)



MSVFNNSEVKYFLLIGIPGLEIAHAWISIPICLMYLIAIVGNCTIIFVIKTEPSLHEPMYYFLTMLAVSD
MGLSFSSLPMTLKI FCFNAMQISPNA CFAQE FFIHVFTAMESSVLLIMSLDRFLATHNPLRYSSILTGSR
VAKIGLILAIIRSTALVLPFPFTRKLKYCKNLLSHLYCLHQDVMKLACSDNKINFIYGFFAALCTMLDF
ALILMSYMLIKLTVLSIASLVERLKALNTCVSHICAVLIFYVPIITLAAIHRFAKHKSPLLVILIAMFL
LVLPMLNPYVICIKTROIREKVLGKLVNLCVR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM 001001279

ORF Size: 936 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001279.1](#), [NP_001001279.1](#)

RefSeq Size: 939 bp

RefSeq ORF: 939 bp


Locus ID: 365324

Cytogenetics: 1q32

MW: 35.3 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RR206275