

Product datasheet for **RR206277L4V**

Fxyd6 (NM_022005) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fxyd6 (NM_022005) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fxyd6
Synonyms:	Php
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022005
ORF Size:	282 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206277).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022005.2 , NP_071288.1
RefSeq Size:	1784 bp
RefSeq ORF:	285 bp
Locus ID:	63847
UniProt ID:	Q91XV6
Cytogenetics:	8q22



[View online »](#)

Gene Summary:

This reference sequence was derived from multiple replicate ESTs and a deposited cDNA, and validated by similar human genomic sequence. This gene encodes a member of a family of small membrane proteins that share a 35-amino acid signature sequence domain, beginning with the sequence PFXVD and containing 7 invariant and 6 highly conserved amino acids. The approved human gene nomenclature for the family is FXVD-domain containing ion transport regulator. FXVD2, also known as the gamma subunit of the Na,K-ATPase, regulates the properties of that enzyme. FXVD1 (phospholemman), FXVD2 (gamma), FXVD3 (MAT-8), FXVD4 (CHIF), and FXVD5 (RIC) have been shown to induce channel activity in experimental expression systems. Transmembrane topology has been established for two family members (FXVD1 and FXVD2), with the N-terminus extracellular and the C-terminus on the cytoplasmic side of the membrane. This gene product, FXVD6, is novel and has not been characterized as a protein. The name "phosphohippolin" has been used in GenBank, but there is no evidence yet of protein phosphorylation. [RefSeq curation by Kathleen J. Sweadner, Ph.D., sweadner@helix.mgh.harvard.edu., Dec 2000]