

Product datasheet for **RR206287**

Spock3 (NM_001107310) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spock3 (NM_001107310) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Spock3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR206287 representing NM_001107310 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTCAAGGTGTCAGCCTTATTGTGTGTGTGTGCAGCGCTTGGTGCAGCCAGACTCTTGCAGCTGCTG
CGGCGGTGGCTGTGGCCGGGGGCGGTTCGGACGGCGGTAAATTTCTGGATGAAAAACAATGGCTCACCAC
AGTCTCTCAGTATGACAAGGAAGTCGGACAGTGAACAATTCGAGACGTAAGTGATGAAGTAGAGGAT
GATTATTTCCGCACTTGAATCCTGGAAAACCAATTGATCAAGCTTTAGATCCTGCTAGGGACCCATGCT
TAAAGACGAAATGCAGCCGCCACAAAGTCTGCATTACTCAAGATGCTCAGACTGCAGTGTGCATCAGTCA
CCGGAGGCTCACACACAGTATGAAAGAAGTAGGAGGAAGCCATAAGCAATGGAGAGGTCTCCCGTCATCT
ACCTGTAAGCAGTGCCCTATAGCCTATGCCAGCCCGGTCTGTGGTCCGATGGGCATTCTTACTCTTCAC
AGTGCAAACTAGAATATCAGGCATGTGTCTTAGGAAAACAGATCTCTATCAAATGTGAAGGGCGATGCC
GTGTCCCTCCGATAAGTCCATGAGCACAGGCAGAAATGTAAGAGAGCCTGCAGTGACCTGGAATTCAGA
GAAGTGGCAAACAGACTGCGAGACTGGTTCAAGGCGCTTCATGAGAGTGAAGCCAGAACAAGAGGACAA
AAGCATTGCTGAGGCCTGAGAGAAGCAGTTTGATACAGTATTTTGCCAATTTGCAAGGATTCACTTGG
CTGGATGTTTAATAGACTTGATACAACTACGACCTGTTATTGGACAGTCAGAGCTTGGCAGCATTTAC
CTTGATAAGAATGAGCATTGCACCAAGGCGTTCTTCAACTCTTGATACATACAAGGACAGTTTGATAT
CTAACAAATGAATGGTGTCTACTGCTTCAGAGACAGCAAGACCCACCTTGCCACACTGAGCTCAGCAAAAT
TCAGAAGCGACAAGGGATAAAGAAGCTCCTAGGACAGTACCTCCCACTCTGTGATGAAGATGGTTTCTAC
AAGCCAACACAATGCCATGGAAGTGTAGGGCAGTGCTGGTGTGTGACAGATATGAAATGAAGTTGTGG
GATCTAGAATAAATGGTGTGTCAGATTGTGCTATAGACTTTGAGATTTCTGGAGACTTTGCAAGTGGAGA
TTTTCTGTAATGGACTGATGATGAAGGTGAAGAAGATGACATGAATGATAAGGATGATATTGAAGACGAT
GATGAAGATGAAGGGATGATGATGATGATGGGATGACCACGATGGATATATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR206287 representing NM_001107310
 Red=Cloning site Green=Tags(s)


MLKVSALLCVCAAWCSQTLAAAAVAVAGGRSDGGNFLDEKQWLTTVSQYDKEVGQWNKFRDVSDEVED
 DYFRTWNPCKPIDQALDPARDPCLKTKCSRHKVCITQDAQTALCISHRRLTHSMKEVGGSHKQWRGLPSS
 TCKQCPIAYASPVCGSDGHSYSSQCKLEYQACVLGKQISIKCEGRCPSPDKSMSTGRNVKRACSDLEFR
 EVANRLRDWFKALHESGSQNKRTKALLRPERSRFDTSILPICKDSLGMFMNRLDTNYDLLLDQSELGSIY
 LDKNEHCTKAFFNSCDTYKDSLISNNEWCYCFQRQDPPCHTELSKIQRQGIKKLLGQYLPLCDEEDGFY
 KPTQCHGSVGCWCVDYRNEVVGSRINGVADCAIDFEISGDFASGDFREWTDDEGEEDDMNDKDDIEDD
 DEDEGDDDDGDDHDGYI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001107310

ORF Size: 1314 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001107310.1, NP_001100780.1

RefSeq Size: 3142 bp


RefSeq ORF: 1317 bp

Locus ID: 306404

Cytogenetics: 16p13

MW: 49.3 kDa

Product images:



Circular map for RR206287