

Product datasheet for **RR206305**

Pwwp2a (NM_001127296) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pwwp2a (NM_001127296) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pwwp2a
Synonyms:	RGD1305891
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RR206305 representing NM_001127296
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGCGGCGCTGGCTGCTGAGGCGGCGGACCGCGCGTCCCCGGGGACGGGGCGCCGGTGAGGCCG
 AGCCGGAGCTAGAGCCATTCCGGGCAGTGAGGCGGCACCCCCCTCCCGGTACAGGCCACGGAGGCGGC
 GGTGCCGACGCGAGGCCGATGGGCGCAGTCCGCCCTCAGGCCGACGAGCAGCCGCTGCCGCCGCCG
 CCGCGCCGCCACCGCCGGGGAGCTGGCTGGTAGCTCCGAGGCTGAGGAGGCGAAGCCGCCTGAGCCCG
 CGGCTGTCCCGTCCGCCGCCGAGCAGCCTCCGGCGCCCGGAGCAGCCGAGGACGCGCCACGGCC
 GCCGCCAGCTCCGGCGCTCGTCCGCCGCGGGCGGGGACTCAGCGGTATCGCACCTCATCCGGGCTCG
 GAGGTGCGGGTACGCTGGACCACATCATCGAGGACGCGCTCGTGGTGTGTTCCGCTCGGGGAGAAGC
 TCTTCTCGGGGCTCATGACCTGTCTAAAAGTTTGGCCCCATGGGATCCCTGTGACAGTATTTCC
 CAAAAGGAATGTAAGGACAAACCGGACCCATGCAGTTCAAAGTAACACATTCCAAGAAGGGATAGAA
 GTCAAGCGGGAAGTGAATGGTGCTATTCCCGATGCCCTTCTCCAGTCCCTCCTCTGATTGCCTGTGGA
 CTTCAAACACCTCCTCTCTTCATGAAGGAGCACCTTATCCTCCCCCTTGTTATCAGGGACACATA
 TAACCAATCAATACCTCAGCCTCCTCCCCGAAAATTAAGCGACCCAAACGAAAATGTACAGGGAAGAA
 CCCACTTCAATAATGAATGCGATAAAGCTACGTCCAGGCAAGTCCTGTGTGATAAGGTAAAAACAGTG
 TTGTTGCAGAAAAGAAGGAAATTAGGAAAGGTAGCAGTGACTCTTCTAGGTATGAAGATAAAAAACGGAG
 AAACGATAGTGTAGTACTGTGAACAAAACTGAAACTGACCATAAAGTGGATGGGAAAAACAAAAC
 GAAAGCCAGAGAAGGAATGCTGTGGTGAAGGTTTCCAGTATCCCTCATAGCAGGGGACAGTAGTCAAGG
 TTTCTGCTCAGGCAAATACATCAAAGCTCAGTTAAATACCAAGAAAGTGCTTCAGAGCAAGAACATGGA
 TCACGCAAAAGCTCGGGAAGTATTGAAAATTGCCAAAGAAAAGGCACAGAAGAAGCAAAGTGAAACCTCC
 ACTTCAAAAATGCGCACTCGAAAGTCCATTTACACGCCGATCAGAATCCTAGCTCAGGTTCTCTCC
 CACCTCGAGTGCCTTAAACCAACAACGGTACAGGAATGAAGAAAACGACTCTTCTCTGAAGACGGGACT
 GGAGAAAATTCGAGTGGCAAGCTGGCCCTAAGCCGAGTCTCGATGTACCTCCACCCGCTCAGCAGGT
 CTCAACAAATGGCAGCTATTACATCAGACAGTGACGAGTCTGCTGCTCCCTTACAGTGTCTGACAGACC
 ACTGTGGATTGAGTGGGAGCATTGAAGTTAACAGTGAAGCGGGCAGCACAAGACAT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR206305 representing NM_001127296
 Red=Cloning site Green=Tags(s)


MAAVAAEAAATAASPGDGGAGEAEPELEPIPGSEAGTLPVPTATEAAVPDGEADGRQSAPQADEQLPPP
 PPPPPGELAGSSEAEAKPPEPAAPVPPPEQPPAAPEQPEDAPRPPAPALVPPAGGDSAVSHLIPGS
 EVRVTLDHIIEDALVVSFRLGEKLFSGVMDLSKRFPGHIGPVTVPKRECKDKPDAMQLQSNTFQEGIE
 VKREVNGAIPDALSPVPPDCLWTSKPPPLFHEGAPYPPPLFIRDYTNQSIPQPPPRKIKRPKRKMYREE
 PTSIMNAIKLRPRQVLCCKNSVVAEKKEIRKGSSDSSRYEDKKRRNDSVATVNKKLTDHKVDGKNQN
 ESQRRNAVVRVSSIPHSRGRVVKVSAQANTSKAQLNTKKVLQSKNMDHAKAREVLKIAKEKAQKKQSETS
 TSKNAHSKVHFTRRYQNPSSGSLPPRVRLKPQRYRNEENDSSLKTGLEKIRSGKLAPKPQSRCTSTRSAG
 LNKWQLLHQVTSPAAPLQCLTDHCGFRLGALKLTVKRAAQRH

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001127296

ORF Size: 1599 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001127296.1](#), [NP_001120768.1](#)

RefSeq Size: 2897 bp


RefSeq ORF: 1602 bp

Locus ID: 303060

Cytogenetics: 10q21

MW: 57.8 kDa

Product images:



Circular map for RR206305