

Product datasheet for **RR206312**

Aga (NM_001031641) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Aga (NM_001031641) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Aga
Synonyms: MGC112709
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR206312 representing NM_001031641
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGCGGAAGTGAATCTCCCTTTCTCTCCTGCCGCTCGTTCTGGGTATACCCCTGGTGCGGGGCT
CCAACCCTCTACCCTGGTCGCAACACTTGGCCTTTAAGAATGCCACCGAAGCAGCGTGGTGGACATT
AGTGTCTGGAGGTTCTGCCCTGGATGCGGTGGAGAAGGGCTGTGCTATGTGTGAAAAGGAGCAATGTGGT
GGTACTGTAGGCTTTGGAGGAAGTCCTGATGAAGTTGGGGAAACCACCCTGGATGCCATGATCATGGATG
GCACTGCTATGGATGTGGGAGCAGTGGGAGGCCTTAGAAGAATTAATAATGCTATTGGAGTGGCTCGGAA
AGTCTTGGAGCATACCACGCACACGCTTTTAGTGGGGACTCAGCCACCAAGTTTGTGTAAGTATGGGG
TTTACCAGTGAGGACTTATCTACCAATACCTCAAGAGCTTTCATTTCAGATTGGCTTTCTCGAAACTGCC
AGCCAAATTATTGGAGAAACGTTATACCAGATCCTTCAAATACTGTGGACCCTACAAACCACCTGATTT
CTTAGAGCAGAATAACCGTGCCCAAAAGAAGTGGATATCCACAGTCATGATACTATTGGCATGGTTGTA
ATCCATAAGACGGGACATACTGCTGCTGGCACATCTACAAATGGTTTAAAATTCAAATACTGGTCGTG
TAGGAGACTACCAATCCCTGGAGCTGGAGCCTATGCTGATGACATGGCTGGAGCAGCTCGGCCACTGG
AGATGGTGACACACTCTGAGATTTCTGCCAAGCTACCAAGCTGTAGAATATATGAGAGGAGGAGATGAT
CCAGCCAGAGCTTGCCAAAAAGTATTTCAAGAATTCAGAAATATTATCCAAAATTTCTTTGGAGCTGTTA
TATGTGCCAATGTGACAGGAAGTATGGTGTGCTGCTGCAACAGACTTCCAACATTTACTCAATTTAGTTT
CATGGTTTATACTCTCTACACAATCAGGCGATTGAAGAAAAAGTAGACTGCATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA




[View online »](#)

Protein Sequence: >RR206312 representing NM_001031641
 Red=Cloning site Green=Tags(s)


MARKWNL P F L L L P L V L G I P L V R G S N P L P L V V N T W P F K N A T E A A W W T L V S G G S A L D A V E K G C A M C E K E Q C G
 G T V G F G G S P D E V G E T T L D A M I M D G T A M D V G A V G G L R R I K N A I G V A R K V L E H T T H T L L V G D S A T K F A V S M G
 F T S E D L S T N T S R A L H S D W L S R N C Q P N Y W R N V I P D P S K Y C G P Y K P P D F L E Q N N R A H K E V D I H S H D T I G M V V
 I H K T G H T A A G T S T N G L K F K I P G R V G D S P I P G A G A Y A D D M A G A A A T G D G D T L L R F L P S Y Q A V E Y M R G G D D
 P A R A C Q K V I S R I Q K Y Y P K F F G A V I C A N V T G S Y G A A C N R L P T F T Q F S F M V Y N S L H N Q A I E E K V D C M

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:
Cloning Scheme:



Plasmid Map:



ACCN: NM_001031641
ORF Size: 1035 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001031641.1 , NP_001026811.1
RefSeq Size:	1215 bp
RefSeq ORF:	1038 bp
Locus ID:	290923
UniProt ID:	P30919
Cytogenetics:	16p11
MW:	37.2 kDa
Gene Summary:	Cleaves the GlcNAc-Asn bond which joins oligosaccharides to the peptide of asparagine-linked glycoproteins.[UniProtKB/Swiss-Prot Function]