

Product datasheet for **RR206336**

Cmtm6 (NM_001007802) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cmtm6 (NM_001007802) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cmtm6
Synonyms: Cklfsf6; Da2-17; LRRGT00102
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR206336 representing NM_001007802
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGAACGGAGCGGTGTACAGTCCCACCACGGAGGCGGCCCGGGCGCTGGCAGGGGCGCTCGCAGCG
GCCTGGCCGCTACTTCTTCTGGGCGGCTCCCATGGTATCGGCGCATCCTCAAGGGCTTCGAGCTGCT
GCTGTCTCTCTGGCCTTCATCTGTGAAGAGGTTGTGTCCGAGTGTGGGCTGTGTGGAGGCTGTACTTC
TTTGAGTTTGTGAGCTGCAGTGCCTTTCTCCTGAGCCTTCTCCTGCTGATCGTGTACTGCACAGCAGTGT
ACGACAGAGTGGATACTGAAAAGTCAAGTCATCGGATTTTACATCACCTGGGCACAGGTGTGTGTT
TCTGCTGGCTTCTATCATCTTCTCTACCCATTCCGGGACCTCGGCTGAGATTGCTGCAATTGTGTTT
GGCTTCTAGGAAGCTCAATGTTTCTACTGGACTTCGTCGTCATGCTGTGTGAGAAGCTGCGGGAGAGCC
CTCTTAGAAAGCCCGAGAACAACGCTAGAGCTGAAGCCCTCACTGAACCACTGAATGCT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR206336 representing NM_001007802
Red=Cloning site Green=Tags(s)

MENGAVYSPTTEAAPGAGRGARSLAAYFFLGRLPWYRRILKGLQLLLSLLAFICEEVVSECGLCGGLYF
FEFVSCSAFLLSLLLLIVYCTAVYDRVDTGKVKSSDFYITLGTGCVFLLASIIIFVSTHSGTSAEIAAIVF
GFLGSSMFLDFVVMLCEKLRRESPLRKPENNRARAEALTEPLNA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_001007802

ORF Size: 549 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007802.3</u> , <u>NP_001007803.1</u>
RefSeq Size:	598 bp
RefSeq ORF:	552 bp
Locus ID:	316035
Cytogenetics:	8q32
MW:	19.9 kDa