

## Product datasheet for **RR206341L4V**

### Dync1i2 (NM\_053880) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dync1i2 (NM_053880) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dync1i2
Synonyms:	Dnci2; Dncic2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053880
ORF Size:	1914 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053880.2</a> , <a href="#">NP_446332.2</a>
RefSeq Size:	2821 bp
RefSeq ORF:	1917 bp
Locus ID:	116659
Cytogenetics:	3q22
Gene Summary:	intermediate chain of cytoplasmic dynein, a retrograde microtubule motor thought to be involved in organelle transport [RGD, Feb 2006]



[View online »](#)