

Product datasheet for **RR206347L3V**

Rgcc (NM_054008) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rgcc (NM_054008) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rgcc
Synonyms:	Rgc32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_054008
ORF Size:	411 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206347).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_054008.1 , NP_446460.1
RefSeq Size:	889 bp
RefSeq ORF:	414 bp
Locus ID:	117183
UniProt ID:	Q9Z2P4
Cytogenetics:	15q12
Gene Summary:	increased expression is detected in oligodendrocytes in response to complement activation; may play a role in cell cycle activation [RGD, Feb 2006]



[View online »](#)