

Product datasheet for **RR206393L4V**

Ebag9 (NM_001009665) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ebag9 (NM_001009665) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ebag9
Synonyms:	MGC94627
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001009665
ORF Size:	639 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206393).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009665.1 , NP_001009665.1
RefSeq Size:	1978 bp
RefSeq ORF:	642 bp
Locus ID:	299864
UniProt ID:	Q5PQP2
Cytogenetics:	7q31
Gene Summary:	May participate in suppression of cell proliferation and induces apoptotic cell death through activation of interleukin-1-beta converting enzyme (ICE)-like proteases.[UniProtKB/Swiss-Prot Function]



[View online »](#)